

Genotyping of rs3824662 polymorphism in GATA3 gene by allele-specific PCR

Nafiul Huda, Md. Ismail Hosen, Pankaj Kumar Sarkar, A.K.M. Mahbub Hasan, A.H.M. Nurun Nabi

Abstract

To amplify the target regions of GATA3 for studying rs3824662 polymorphism, a web-based primer designing tool available at <http://bioinfo.ut.ee/primer3-0.4.0> facilitates crafting the set of outer primers [forward outer primer: 5'-TTGCAAATGGAAGAGGGTCT-3' and reverse outer primer: 5'-ACCCTGCAAATGAGAGGAAA-3'] and inner primers [G specific primer: 5'-TGAGATTAAACACAAACACGtTG-3' and T specific primer: 5'-CTGAGATTAAACACAAACACGaTT-3'] to perform allele-specific PCR that specifically amplify the GATA3 gene. G/T allele-specific primer is very sensitive and binds only with the DNA fragments containing G/T nucleotide template DNA. Forward outer and reverse outer primers will amplify a specific 691 base pair region of GATA3. This 691 base pair region contains the desired polymorphic site of interest. Allele-specific primers produce 506 bp products upon the presence of either a G allele or T allele. PCR reaction in a total volume of 15 µL with an initial denaturing step of 5 min at 95°C, followed by 40 cycles of 30 s at 95°C, 45 s at 55°C, and 1 min at 72°C, and a final extension step of 5 min at 72°C will produce better reproducible result. The reaction mixture contains both forward outer and reverse outer primers and contains either the G allele-specific forward inner primer or the T allele-specific forward inner primer. Primer concentration used was 200 nM each. When the G allele is present in a specific DNA sequence, only the 506 base pair band will be found upon gel electrophoresis of PCR amplicons with the G allele-specific primer. In this case, no such 506 base pair band will be obtained upon gel electrophoresis with the T allele-specific primer. On the other hand, for mutant TT genotypes, the opposite phenomenon will occur when the T allele is present, which will be reflected by the presence of the 506 base pair DNA band. In the case of the heterozygous genotype (GT), both inner primers for G and T alleles will bind with specific DNA sequences to produce 506 base pair bands.


Citation: Nafiul Huda, Md. Ismail Hosen, Pankaj Kumar Sarkar, A.K.M. Mahbub Hasan, A.H.M. Nurun Nabi Genotyping of rs3824662 polymorphism in GATA3 gene by allele-specific PCR. **protocols.io**
dx.doi.org/10.17504/protocols.io.k45cyy6


Published: 04 Dec 2017


Before start

Take precaution as required for a general PCR reaction.

Materials

 GoTaq(R) G2 Hot Start Green Master Mix,
1,000 reactions M7423 by Promega

 Agarose, LE, Analytical Grade, 500gm v3125
by Promega

 Ethidium bromide [EB, EtBr] EB0195.SIZE.25g
by Bio Basic Inc.

✓ primers by Contributed by users

Protocol

Primer design

Step 1.

The set of outer primers [forward outer primer: 5'-TTGCAAATGGAAGAGGGTCT-3' and reverse outer primer: 5'-ACCCTGCAAATGAGAGGAAA-3'] and inner primers [G specific primer: 5'-TGAGATTAAACACAAACACGtTG-3' and T specific primer: 5'-CTGAGATTAAACACAAACACGaTT-3'] were designed.

PCR condition

Step 2.

G/T allele-specific primer is very sensitive and binds only with the DNA fragments containing G/T nucleotide template DNA. Forward outer and reverse outer primers will amplify a specific 691 base pair region of GATA3. This 691 base pair region contains the desired polymorphic site of interest. Allele-specific primers produce 506 bp products upon the presence of either a G allele or T allele.


PCR reaction in a total volume of 15 µL with an initial denaturing step of 5 min at 95°C, followed by 40 cycles of 30 s at 95°C, 45 s at 55°C, and 1 min at 72°C, and a final extension step of 5 min at 72°C will produce better reproducible result. The reaction mixture contains both forward outer and reverse outer primers and contains either the G allele-specific forward inner primer or the T allele-specific forward inner primer. Primer concentration used was 200 nM each.

Output

Step 3.

When the G allele is present in a specific DNA sequence, only the 506 base pair band will be found upon gel electrophoresis of PCR amplicons with the G allele-specific primer. In this case, no such 506 base pair band will be obtained upon gel electrophoresis with the T allele-specific primer. On the other hand, for mutant TT genotypes, the opposite phenomenon will occur when the T allele is present, which will be reflected by the presence of the 506 base pair DNA band. In the case of the heterozygous genotype (GT), both inner primers for G and T alleles will bind with specific DNA sequences to produce 506 base pair bands.

 EXPECTED RESULTS



Warnings

1. As T_m value is critical so to run the PCR reaction, one should set its own T_m value considering ± 2 degree centigrade temperature.
2. Please use gloves when someone is using Ethidium bromide as it is a carcinogenic compound. Also, disposal of this compound should be taken care of.