



24 Hour Metabolic Response to LPS in Monocytes

Brandt Pence¹

¹University of Memphis

dx.doi.org/10.17504/protocols.io.ufdeti6

Oct 08, 2018

Working



Brandt Pence

University of Memphis




PROTOCOL STATUS

Working

We use this protocol in our group and it is working

MATERIALS

NAME	CATALOG #	VENDOR
 Cell-Tak	354240	Corning
 Seahorse XFp FluxPak	103022-100	Agilent Technologies
 Seahorse Base Medium DMEM	102353-100	Agilent Technologies
 100 mM Sodium Pyruvate	S8636	Sigma
 200 mM L-Glutamine	G7513	Sigma
 45% D-()-Glucose	G8769	Sigma
 0.1 M Sodium Hydroxide Solution	71395	Sigma
 Lipopolysaccharide (1 microgram/ml in LPS stored -20decC)	L4391	Sigma Aldrich
 RPMI-1640 without phenol red	11835055	Gibco - Thermo Fisher
 Fetal Bovine Serum	10437028	Gibco - Thermo Fisher
 Penicillin-Streptomycin	15140122	Gibco - Thermo Fisher
 0.1M Sodium Bicarbonate Buffer pH 8.0 Sterile	View	Contributed by users
 Pipettes and P1000 P200 P10 tips	View	Contributed by users

Day 1 Materials to Prepare

1 Stimulation media

Add 500 µl fetal bovine serum and 50 µl pen-strep to 4.45 ml RPMI, sterile filter

Cell-Tak coated XFp plate.

Prepare Cell-Tak. Add Cell-Tak and sodium hydroxide to 0.1 M sodium bicarbonate so that each well will receive 0.56 µg Cell-Tak, and sodium hydroxide concentration will be 0.63 mM. Cell-Tak concentration varies by batch, so calculations for each new batch will need to be performed. Add 25 µl Cell-Tak solution to each well. Plate can be stored at room temperature (at least 20 minutes) until use (for a few hours). Plates can also be prepped prior to assay day by incubating plates in Cell-Tak (at least 20 minutes), then aspirating remaining solution, air drying, and storing at 4°C until use.

Monocytes

Monocytes should be prepared as directed in the monocyte isolation protocol. Isolated monocytes should be diluted to a concentration of 3×10^6 cells / ml in prepared stimulation media prior to use in the assay. Use monocytes immediately.

Day 1 Procedure

- 2 Aspirate Cell-Tak solution if not already done.
- 3 Add 50 µl medium to wells A and H and 50 µl cells to wells B-G. Samples are generally run in duplicate or triplicate on each plate. Cell number is 1.5×10^5 cells per well.
- 4 Place plate in the carrier and place in centrifuge. Spin $300 \times g$ for 1 minute **without brake**.
- 5 Add 130 µl stimulation medium to each well A-H (final volume 180 µl).
- 6 Prepare LPS: 3 µl of 1 µg/ml LPS in 297 µl stimulation medium
- 7 Add 20 µl media to wells A-D, H and 20 µl LPS dilution to wells E-G
- 8 Incubate plate at 37°C/5%CO₂ incubator for 23 hours.
- 9 Turn on XFp analyzer to warm up overnight
- 10 Hydrate XFp sensor cartridge
Add 200 µl XF calibration solution (included with FluxPak) to each well.
Add 400 µl sterile PBS or H₂O to each moat.
Incubate overnight at 37°C in non-CO₂ incubator.

Day 2 Materials to Prepare

- 11 Seahorse DMEM Media
Add 50 µl pyruvate, 50 µl L-glutamine, 20 µl glucose to 5 ml Seahorse Base Medium. Sterile filter before use.

Day 2 Procedure

- 12 Place plate in the carrier and place in centrifuge. Spin $300 \times g$ for 1 minute **without brake**.
- 13 Aspirate supernatant and replace with 200 µl Seahorse DMEM medium

- 14 Incubate plate 1 hr at 37°C in non-CO₂ incubator
- 15 While plate is incubating, complete steps 16-17.
- 16 Remove sensor cartridge from incubator and remove and reinsert sensors briefly to clear air bubbles.
- 17 Set up program on Seahorse XFp and calibrate sensor cartridge (**remove lid**).
Basal: 6 measurements
- 18 After 1 hour cell incubation, remove utility plate from XFp and insert cell plate (**remove lid**).
- 19 Run Assay.
- 20 After run is completed, image each well by photomicroscopy or collect and isolate protein from each well to normalize cell numbers.

Data Analysis

- 21 Data are analyzed by mean ECAR.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited