OPEN ACCESS



# Quick Protocol for Oligonucleotide Cleanup Using the Monarch® PCR & DNA Cleanup Kit (5 μg) (NEB #T1030)

#### **New England Biolabs**

#### **Abstract**

Quick Protocol for Oligonucleotide Cleanup Using the Monarch® PCR & DNA Cleanup Kit (5 μg) (NEB #T1030)

**Citation:** New England Biolabs Quick Protocol for Oligonucleotide Cleanup Using the Monarch® PCR & DNA Cleanup Kit (5 μg) (NEB #T1030). **protocols.io** 

dx.doi.org/10.17504/protocols.io.nyfdftn

Published: 06 Apr 2018

#### **Guidelines**

**DNA Cleanup and Concentration:** for the purification of up to 5  $\mu$ g of DNA (ssDNA > 200 nt and dsDNA > 50 bp) from PCR and other enzymatic reactions. A <u>detailed protocol</u> and a <u>quick protocol</u> are available for your convenience.

**Oligonucleotide Cleanup (steps)**: for the purification of up to 5  $\mu$ g of DNA fragments  $\geq$  15 bp (dsDNA) or  $\geq$  18 nt (ssDNA). Expected recovery is > 70%. When purifying ssDNA of any size, recovery can be increased by using this protocol; however, it is important to note that this protocol shifts the cutoff for smaller fragments to 18 nt (rather than 50 nt for the DNA Cleanup and Concentration Protocol).

#### Before start

- All centrifugation steps should be carried out at 16,000 x q (~13,000 RPM).
- Add 4 volumes of ethanol (≥ 95%) to one volume of DNA Wash Buffer.

#### **Materials**

Monarch® PCR & DNA Cleanup Kit (5 μg) T1030 by New England Biolabs

#### **Protocol**

#### Step 1.

Add 100 µl DNA Cleanup Binding Buffer to the 50 µl sample.

#### **■** AMOUNT

100 µl Additional info: DNA Cleanup Binding Buffer

#### NOTES

Danielle Freedman 21 Mar 2018

We recommend a sample volume of 50  $\mu$ l. For smaller samples, adjust the volume with nuclease-free water.

#### Step 2.

Add 300  $\mu$ l ethanol ( $\geq$  95%). Mix well by pipetting up and down or flicking the tube. Do not vortex.

## **■** AMOUNT

300 µl Additional info: Ethanol (≥ 95%)

#### Step 3.

Insert column into collection tube and load sample onto column.

#### Step 4.

Spin for 1 minute, then discard flow-through.

#### Step 5.

Re-insert column into collection tube.

#### Step 6.

Add 500 µl DNA Wash Buffer and spin for 1 minute.

#### **■** AMOUNT

500 μl Additional info: DNA Wash Buffer

#### Step 7.

Discard flow-through.

#### Step 8.

**Repeat steps 5-7 (Optional).** Recommended for removal of enzymes that may interfere with downstream applications (e.g., Proteinase K).

# **P**GOTO

Repeating steps 5-7 -> go to step #5

#### Step 9.

**Transfer column to a clean 1.5 ml microfuge tube.** Use care to ensure that the tip of the column does not come into contact with the flow-through. If in doubt, re-spin for 1 minute.

## Step 10.

Add  $\geq$  6  $\mu$ l of DNA Elution Buffer to the center of the matrix.

# **■** AMOUNT

6 μl Additional info: DNA Elution Buffer

# **Step 11.**

Wait for 1 minute, then spin for 1 minute to elute DNA.

#### **P** NOTES

#### Danielle Freedman 21 Mar 2018

Typical elution volumes are  $6\text{--}20~\mu l$ . Nuclease-free water (pH 7-8.5) can also be used to elute the DNA. Yield may slightly increase if a larger volume of DNA Elution Buffer is used, but the DNA will be less concentrated.

# **Warnings**

Please refer to SDS for safety warnings.