

# RNA extraction from Escherichia coli version 3

Alice Pawlowski, Lutz Berwanger

## Abstract

RNA extraction from E. coli cells based on the method described by Chomczynski and Sacchi, 1987

**Citation:** Alice Pawlowski, Lutz Berwanger RNA extraction from Escherichia coli. **protocols.io**

dx.doi.org/10.17504/protocols.io.kzycx7w

**Published:** 27 Nov 2017

## Guidelines

RNA is sensitive to degradation! Wear gloves, keep samples on ice when possible, use filter-tips and RNase free reagents. Pre-cool centrifuges and store isolated RNA-samples immediately at -20 or -80°C.

## Before start

Always keep your samples on ice!

## Materials

Roti®-C/I by [Carl Roth](#)

Roti®-P/C/I by [Carl Roth](#)

## Protocol

### Cell preparation

#### Step 1.

- mix 1 ml of cells with 200 µl 'stopmix'- solution (5 % phenol in ethanol) in a 2 ml tube → stops RNA production in the cells

### SAFETY INFORMATION

#### Wear safety gear

### Cell preparation

## Step 2.

- centrifuge for 5 min at 4°C and 14000 x g
- discard the supernatant and resuspend the pellet in 1 ml NucleoZOL (Macherey and Nagel), place on dry ice
- proceed to next step or store cells at -20 or -80 °C

### TEMPERATURE

4 °C Additional info:

RNA-isolation

## Step 3.

- incubate the sample at 65 °C and 250 rpm (Thermomixer) for 10 min
- mix with 400 µl Chloroform/Isoamylalcohol (Roti®-C/I) by inverting for 10 s

### TEMPERATURE

65 °C Additional info:

RNA-isolation

## Step 4.

- centrifuge at 4°C for 10 min at 14000 x g
- transfer aqueous phase to a new reaction tube, work on ice
- mix with 450 µl Penol/Chloroform/Isoamylalcohol (Roti®-Aqua-P/C/I)
- centrifuge at 4°C for 10 min at 14000 x g
- transfer aqueous phase to a new reaction tube and add 1 Vol. icecold Isopropanol + 20 µl 3 M Na-Acetate (pH 5.2) and mix

### TEMPERATURE

4 °C Additional info:

### SAFETY INFORMATION

**Wear safety gear!**

RNA-isolation

## Step 5.

- leave RNA at least 30 min at -20 °C or store over night

RNA-isolation

## Step 6.

- centrifuge at 4°C for 30 min at 14000 x g
- remove the supernatant (take care of the RNA-pellet) and add 350 µl of icecold 75% ethanol

### TEMPERATURE

4 °C Additional info:

RNA-isolation

## Step 7.

- centrifuge at 4°C for 5 min at 14000 x g

🌡 **TEMPERATURE**

4 °C Additional info:

RNA-isolation

**Step 8.**

- add 350 µl of icecold 75% ethanol
- centrifuge for 5 min at 4°C and 14000g

🌡 **TEMPERATURE**

4 °C Additional info:

RNA-isolation

**Step 9.**

- remove the supernatant and dry the pellet at room temperature for ca.15 min
- resuspend the pellet in 30 µl molecular biology grade water or TE-buffer (10 mM Tris/HCL pH 8.0 , 1 mM EDTA)

## Warnings

Phenol is toxic! Work under the hood, always wear protective gear and change contaminated gear immediately. Collect solid and liquid waste in special waste containers.