# Natural Transformation of Campylobacter jejuni Version 2

#### Jessica Sacher

# **Abstract**

Citation: Jessica Sacher Natural Transformation of Campylobacter jejuni. protocols.io

dx.doi.org/10.17504/protocols.io.magc2bw

Published: 18 Dec 2017

## **Protocol**

## Day 1

## Step 1.

Thaw cells, grow 24-48 hours.

# Day 2

#### Step 2.

Restreak, grow 18-24 h.

#### Day 3

# Step 3.

Early in the day: Streak cells in a thin layer onto BHI + 2% yeast (BHI+Y) (yeast is optional but is helpful). Grow 6 h.

# Day 3

# Step 4.

Pre-warm 1-2 BHI+Y plates.

#### Day 3

#### Step 5.

Harvest 6-h plate of cells in 2 mL PBS, wash 1x in PBS, resuspend in 250 uL PBS.

#### Day 3

# Step 6.

Spot entirety of cell suspension onto pre-warmed BHI+Y plates by filling a pipette tip and spotting discrete 20-uL drops across the plate(s).

#### Day 3

# Step 7.

Pipette plasmid DNA (20 ng/uL in water) atop each spot (10-20 uL drops). Let sit to dry a few mins. Incubate O/N (right-side-up).

#### Day 4

# Step 8.

Pre-warm antibiotic-containing BHI plates (will need 6)

## Day 4

#### Step 9.

Streak growth from cell+DNA spots onto BHI + antibiotic plates in **a thin layer** (streak out as much of the growth as you can - use several plates). Incubate 2-5 days, checking each day for colonies.

# Day 6

# Step 10.

Patch colonies onto a new selective plate and grow overnight.

# Day 7

# **Step 11.**

Perform colony PCR using gene-specific or vector-specific primers to check for insert DNA.

# Day 7

# **Step 12.**

Streak out full plates of growth from at least 3 different positive colonies in order to make frozen stocks.

# NOTES

# Jessica Sacher 16 Dec 2017

Tip: It can be easier to transform 81-176 than 11168, so if you're having trouble getting DNA into 11168, try putting pDNA into 81-176 first, extracting pDNA from those cells, and transforming that into 11168.