

# TBE buffer (1L of 10X)

Bao Thai

## Abstract

**Citation:** Bao Thai TBE buffer (1L of 10X). **protocols.io**

[dx.doi.org/10.17504/protocols.io.mpd5i6](https://dx.doi.org/10.17504/protocols.io.mpd5i6)

**Published:** 16 Jan 2018

## Protocol

### Step 1.

Tris base (121.14 g/mol)	121.1 g	1 M
EDTA (disodium salt dihydrate) mw 372.24 g/mol	7.4 g	0.02 M
Boric acid	61.8 g	1 M
If using for RNA work, prepare with autoclaved water. Store up to 6 months at RT.		