

Skeletonema DNA extraction by Plant DNAzol

Adrian Clarke

Abstract

Citation: Adrian Clarke Skeletonema DNA extraction by Plant DNAzol. **protocols.io**

dx.doi.org/10.17504/protocols.io.g6bbzan

Published: 09 Feb 2017

Protocol

Preparations

Step 1.

- Prepare **Plant DNAzol** by adding RNase (100ug/mL, make 0,6mL/sample)
- Prepare **Wash solution**, 1 part **Plant DNAzol** to 0,75 parts 99% ethanol (2x0,6mL/sample)
- Prepare **75% ethanol** solution (5*0,6mL/sample)

Pellet Skeletonema culture (1500 x g, 5 min), remove growth media and dissolve in f/2 or growth media ~250uL.

Step 2.

Use Swing-out rotor

Put liquid N2 in mortar, drop Skeletonema suspension into N2 and grind once its evaporated.

Step 3.

Transfer ground cells to 2mL eppendorf tube and add 600uL prepared DNAzol (with RNase).

Step 4.

Mix by inversion; incubate at 25°C (RT) for 5 min with a little shake.

Step 5.

Add 600uL Chloroform, mix by vortexing in short bursts, and incubate at 25°C (RT) for 5 min with a little shake.

Step 6.

Centrifuge 10 min at 12000g

Step 7.

Transfer the supernatant (aqueous phase) to a fresh eppendorf tube. (Around 800 uL)

Step 8.

Add 1200uL 99% ethanol to aqueous phase, mix by inversion and incubate for 5 min at RT.

Step 9.

Adjust volume to match 75% of supernatant if using other than 800uL

Pellet DNA by centrifugation, 4 min, 7000g. Remove the supernatant.

Step 10.

Add 600uL Wash solution to precipitate, mix by pipetting,

Step 11.

Incubate samples 5 min at RT.

Step 12.

Pellet DNA by centrifugation, 4 min, 7000g. Remove the supernatant by pipette

Step 13.

Repeat the Washing Procedure one more time (Step 11-

Step 14.

. Add 600uL 75% ethanol solution

Step 15.

Pellet DNA by centrifugation, 4 min, 7000g. Remove the supernatant as much as possible.

Step 16.

Repeat step 15+16 three-four additional times.

Step 17.

Evaporate the remaining ethanol in speed

Step 18.

. Dissolve the DNA in 80uL Water

Step 19.

Centrifuge 12000g for 5 min

Step 20.

Take the supernatant (75uL); add 0,75 ul 10xTE.

Step 21.

Makes for 0.1xTE for storage in 4C or -20