



Sep 25, 2019

## Preparation of LB Media V.2

NUS iGEM<sup>1</sup><sup>1</sup>National University of Singapore

1

Works for me

dx.doi.org/10.17504/protocols.io.7pdhmi6

NUS iGEM  
National University of Singapore

## MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

Water refers to sterilized deionized water

Luria Broth Base (Miller's LB Broth Base)&amp;trade;, powder

12795027

Thermo Fisher

## STEPS MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

Luria Broth Base (Miller's LB Broth Base)&amp;trade;, powder

12795027

Thermo Fisher


Water refers to sterilized deionized water

- 1 Weigh  25 g of Luria Broth Base powder.

Luria Broth Base (Miller's LB Broth  
Base)&trade;, powder

by Thermo Fisher

Catalog #: 12795027

- 2 Add the powder into  1 L of water.



Water refers to sterilized deionized water

- 3 Autoclave entire bottle of LB media.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited