



Oct 18, 2019

E. coli and B. subtilis Colony PCR

[iGEM Dusseldorf¹](#)¹Heinrich-Heine Universität Düsseldorf[1](#) *Works for me* [dx.doi.org/10.17504/protocols.io.8fhhtj6](https://doi.org/10.17504/protocols.io.8fhhtj6)

iGEM Dusseldorf

- 1 Prepare a PCR master mix for 9 reactions, as follows:

	1x Reaction (Volume; μ L)	9x Reactions (Volume; μ L)
Sterile MilliQ Water	8	72
Red Taq 2x Master Mix(1.5 mM MgCl ₂)	10	90
F Primer (10 μ M)	1	9
R Primer (10 μ M)	1	9
Total Volume	20	180

- 2 Aliquot 20 μ L of the master mix in 8 tubes of a PCR strip.
- 3 Transfer cells from a single colony (8 in total) using a sterile P2 pipette tip or inoculating loop

- 4 Run thermocycler using the following program (extension for 1.5 kb insert):

Step	Temp (°C)	Time (mm:ss)	Purpose
1	95	3:00	Initial Denaturation
25 cycles of Steps 2 to 4:			
2	95	0:20	Denaturation
3	57	0:30	Annealing
4	72	1:00	Extension
5	72	5:00	Final Extension
6	4	∞	Storage



Extension time may vary depending on insert size

- 5 Run samples on a gel at 80-150 V until the dye line is approximately 75-80% of the way down the gel. A typical run time is about 1-1.5 hours, depending on the gel concentration and voltage.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited