

Dec 31, 2019

Cloning shRNA Oligos into pLKO.1 V.2

Addgene The Nonprofit Plasmid Repository¹¹Addgene

1

Works for me

dx.doi.org/10.17504/protocols.io.bawmifc6



Addgene The Nonprofit Plasmid Repository



ABSTRACT

This is the protocol accompanying the "pLKO.1 – TRC Cloning Vector". For information about the pLKO.1-TRC cloning vector and tips on designing shRNA oligos for pLKO.1 see Addgene's website: <http://www.addgene.org/tools/protocols/plko/>

EXTERNAL LINK

<http://www.addgene.org/tools/protocols/plko/>

STEPS MATERIALS

NAME	CATALOG #	VENDOR
NEBuffer 3 - 5.0 ml	B7003S	New England Biolabs
NEBuffer 1 - 5.0 ml	B7001S	New England Biolabs
Agel - 300 units	R0552S	New England Biolabs
EcoRI - 10,000 units	R0101S	New England Biolabs
EcoRI - 10,000 units	R0101S	New England Biolabs
T4 DNA Ligase Reaction Buffer - 6.0 ml	B0202S	New England Biolabs
T4 DNA Ligase - 20,000 units	M0202S	New England Biolabs

Annealing Oligos





- 1 Resuspend oligos in ddH₂O to a concentration of 20 µM.
- 2 Add 5µl Forward oligo
5 µl
- 3 Add 5µl Reverse oligo
5 µl
- 4 Add 5 µL 10x NEB buffer 2
5 µl



NEBuffer 3 - 5.0 ml

by New England Biolabs

Catalog #: B7003S

- 5 Add 35 μL ddH₂O
 **35 μl**
- 6 Incubate for 4 minutes at 95°C in a PCR machine or in a beaker of boiling water.
 **00:04:00**
- 7 Incubate the sample at 70°C for 10 minutes in a PCR machine.
 **00:10:00**
- 8 Slowly cool to room temperature over the period of several hours.
 **03:00:00**





This will take a few hours, but it is important for the cooling to occur slowly for the oligos to anneal.




If using a beaker of water, remove the beaker from the flame, and allow the water to cool to room temperature.

Digesting pLKO.1 TRC Cloning Vector

- 9 Mix: 6 μg pLKO.1 TRC-cloning vector (maxiprep or miniprep DNA)
 **6 μg**
- 10 with 5 μL 10x NEB buffer 1
 **5 μl**





NEBuffer 1 - 5.0 ml
by New England Biolabs
Catalog #: B7001S

- 11 with 1 μL Agel
 **1 μl**



Agel - 300 units
by New England Biolabs
Catalog #: R0552S

- 12 bring to 50 μL ddH₂O
 **50 μl**
- 13 Incubate at 37°C for 2 hours.
 **02:00:00**
- 14 Purify with Qiaquick gel extraction kit, eluting in 30 μL of ddH₂O.
- 15 Digest eluate with EcoRI by mixing: 30 μL pLKO.1 TRC-cloning vector digested with Agel

16 with 5 μ L 10x NEB buffer for EcoRI

 5 μ L



EcoRI - 10,000 units

by New England Biolabs

Catalog #: R0101S

17 with 1 μ L EcoRI

 1 μ L



EcoRI - 10,000 units

by New England Biolabs

Catalog #: R0101S

18 and 14 μ L ddH₂O

 14 μ L

19 Incubate at 37°C for 2 hours.

 02:00:00

20 Run digested DNA on 0.8% low melting point agarose gel until you can distinctly see 2 bands, one 7kb and one 1.9kb.



When visualizing DNA fragments to be used for ligation, use only long-wavelength UV light. Short wavelength UV light will increase the chance of damaging the DNA.

21 Cut out the 7kb band and place in a sterile microcentrifuge tube.

22 Purify the DNA using a Qiaquick gel extraction kit. Elute in 30 μ L of ddH₂O.

23 Measure the DNA concentration.

Ligating and Transforming into Bacteria

24 Use your ligation method of choice. For a standard T4 ligation, mix: 2 μ L annealed oligo from "Annealing Oligos" section above.

 2 μ L

25 With 20 ng digested pLKO.1 TRC-cloning vector from the "Digesting pLKO.1 TRC Cloning Vector" section above.

 20 ng



If you were unable to measure the DNA concentration, use 1 μ L

26 With 2 μ L 10x NEB T4 DNA ligase buffer

 2 μ L



T4 DNA Ligase Reaction Buffer - 6.0 ml

by New England Biolabs

Catalog #: B0202S

27 With 1 μ L NEB T4 DNA ligase

 1 μ L



T4 DNA Ligase - 20,000 units

by New England Biolabs

Catalog #: M0202S

28 Bring up to 20 μ L with ddH₂O

29 Incubate at 16°C for 4-20 hours.

 04:00:00

30 Transform 2 μ L of ligation mix into 25 μ L competent DH5 alpha cells, following manufacturer's protocol.

31 Plate on LB agar plates containing 100 μ g/mL ampicillin or carbenicillin (an ampicillin analog).



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited