



Jan 06,
2020

El-cheep-o miniprep

Forked from [El-cheep-o bacterial DNA isolation / miniprep](#)

[Magdalena Julkowska](#)¹

¹King Abdullah University of Science and Technology

1 *Works for me* dx.doi.org/10.17504/protocols.io.peudjew

Salt Lab KAUST



Magdalena Julkowska
King Abdullah University of Science and Technology

ABSTRACT

This is a protocol for cheep miniprep isolation. Be aware that the DNA quality might not be sufficient to use this protocol for cloning purposes. It works very well for the digestion analysis though.

BEFORE STARTING

Prepare following buffers and keep P3 at 4C:

P1: into 250 ml H₂O

EDTA 0.093g 1mM

TRIS 1.515g 50mM(pH 8.0)

AUTOCLAVE

P2: into 250 ml H₂O

NaOH 2g 0.2N

SDS 2.5g 1%

AUTOCLAVE

P3: into 250 ml H₂O

K-acetate 73.63g

Acetic Acid 27.5ml

pH 5.5

AUTOCLAVE

- 1 · Spin 3ml of the overnight culture to pellet cells (1.5-ml for time) @13000 rpm for 5 min
- 2 · Remove the supernatant
- 3 · Add 250 ml of sol P1+Rnase 5ml and vortex to resuspend cells
- 4 · Take 15 min @ room temperature
- 5 · Add 250 ml of sol P2 and mix gently
- 6 · Take 5 min @ room temperature

- 7 ·Add 350 ml of sol P3 and mix gently
- 8 ·Leave 10 min on ice
- 9 ·Spin @ the max speed for 10 min
- 10 ·Take 750 ml of the supernatant and put in a new eppendorf tube and add in this tube 750 ml of isopropanol
- 11 ·Mix well
- 12 ·Spin for 30 min @ 14000 rpm
- 13 ·Remove the supernatant
- 14 ·Spin again for other 10 min
- 15 ·Take @ 37°C for 10 min
- 16 ·Add 50 ml of H₂O
- 17 ·Store @ -20°C



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited