

# Poultry Enterprise Environmental Sample Handling and Processing for Salmonella Spp. Detection

Helen Crabb

## Abstract

This protocol details the methodology for processing environmental samples for Salmonella spp. detection.

**Citation:** Helen Crabb Poultry Enterprise Environmental Sample Handling and Processing for Salmonella Spp. Detection. protocols.io

dx.doi.org/10.17504/protocols.io.n6pdhdn

**Published:** 01 Apr 2018

## Materials



Buffered Peptone Water CM0509 by Oxoid  
Microbiology Products - Thermo Fischer



Nasco Whirl Pak B00736 by Contributed by  
users



Nasco Whirl-Pak B01297 by Contributed by  
users



Mirella Research Sterile Bags (400 x 500  
mm) AAF2 by Contributed by users

## Protocol

### Primary sample collection and handling

#### Step 1.

##### Primary Sample Collection

The sampling methodology is described in the protocol Poultry Enterprise Environmental Sampling Methodology 10.17504/protocols.io.n6mdhc6.

##### Primary Sample Handling

All samples are collected into a sterile plastic container or bag for processing (Whirl-pak)

All primary samples are to be labelled with the date of collection, sample type, sample ID, sample location, and flock ID

All wet primary samples (rinsate or swabs pre-moistened with buffered peptone water) are refrigerated as soon as practical from sample collection and during transportation

All dry collected primary samples (chick papers, hatch debris) are stored cool and out of the sun prior to and during transportation

## Sample Processing

All samples were delivered to the laboratory and processed on the same day as collection.

### Primary Sample Processing

To maintain sample integrity, or minimise cross contamination during processing

All samples collected into a sterile plastic container for processing were kept in the original collection container or bag for processing, unless described otherwise.

The primary sample is processed as described in the protocol below

### Microbiological testing of the primary sample

Microbiological testing of the primary sample is detailed in the protocol Poultry Enterprise Environmental Microbiological Testing for Salmonella spp. [10.17504/protocols.io.n6sdhee](https://10.17504/protocols.io.n6sdhee)

## Primary Samples

### Step 2.

The primary sample types collected from each location, material used for sample collection and the sample pool sizes for microbiological testing are summarised in the table below.

Location	Sample Type	Sample Material	Primary Sample Pool Size (No. Units or Volume)
Parent Egg Production	Dust	Gauze Swab	4
	Manure Belt	Gauze Swab	4
	Egg Belt	Gauze Swab	4
	Boot Swab	Boot Swab	1
Parent Rearing	Drag Swab	Swab Assembly	1
	Chick Papers	Paper or Cardboard	10
Hatchery	Chick Papers	Paper	10
	Hatch Debris	Hatch remains: shell and remnants	3 x 100 g (sub-sample of sample)
Broiler Farm	Boot Swab	Boot Swab	1
	Dust - Wall	Gauze Swab	4
	Dust - Fan	Gauze Swab	4-6
Processing	Carcass Rinse	Rinsate	10 mL (sub-sample of rinsate)
	Portion Rinse	Rinsate	10 mL (sub-sample of rinsate)

## Buffered Peptone Water

### Step 3.

Sterile buffered peptone water (BPW) is added as part of the primary sample processing

BPW is mixed accordance with the manufacturers instructions and sterilised by autoclaving

## Chick Papers

### Step 4.

Ten chick papers collected from day old chick deliveries or at the hatchery is pooled as a single sample

900mL of BPW is added to each sample of 10 chick papers

Leave the sample at room temperature to incubate for approximately 30 minutes.

The sample is then manually macerated by massaging the paper within the sterile bag.

Massage gently to ensure the paper is fully wet through and BPW covers the whole sample

**Do not** overmix

The whole sample is incubated

 **AMOUNT**

900 ml Additional info:

Buffered Peptone Water

### Hatch Debris

#### **Step 5.**

Each collected sample of hatch debris is crushed manually and mixed thoroughly by shaking/tipping.

Three 100g sub-samples are aliquoted from each hatch debris sample as the primary sample.

Each 100g sub-sample is processed as a separate sample

100mL of BPW is added to each 100g sample and gently mixed

**Do not** overmix

The whole sample is incubated

 **AMOUNT**

100 ml Additional info:

Buffered Peptone Water

 **AMOUNT**

100 g Additional info:

Hatch Debris

### Carcass or portion rinse samples

#### **Step 6.**

Carcass rinse samples and portion rinse samples are collected in accordance with the Australian Standard AS4465.

#### **Carcass Rinse Sample**

A fresh chicken carcass is collected from the processing line or chiller into a sterile plastic bag.

500 mL of BPW is added to a single bagged chicken, and massaged gently for 2 - 3 minutes.

A sub-sample (10 mL) of the rinsate is collected as the primary sample

## Portion rinse samples

A fresh chicken portion from the processing line into a sterile plastic bag

200 mL of BPW is added to a single portion and massaged gently for 2 - 3 minutes.

A sub-sample (10 mL) of the rinsate is collected as the primary sample

The whole rinsate sample is incubated

### AMOUNT

500 ml Additional info:

Carcass Rinse -Buffered

Peptone Water

### AMOUNT

200 ml Additional info:

Portion Rinse - Buffered

Peptone Water

## Boot Swabs

### Step 7.

Each boot swab is processed as a single sample.

200mL of BPW is added to each boot swab

Massage sample gently to ensure the swabs are fully wet or covered with BPW

**Do not** overmix

The whole sample is incubated

### AMOUNT

200 ml Additional info:

Buffered Peptone Water

## Dust, Manure Belt, Egg Belt , and drag swab assemblies

### Step 8.

Each pool of gauze swabs (4 or 6 swabs or a single drag swab assembly) is processed as a single sample

100 mL of BPW is added to each pool of gauze swabs

Massage gently to ensure the swabs are fully wet and covered by the BPW

**Do not** overmix

The whole sample is incubated

### AMOUNT

100 ml Additional info:  
Buffered Peptone Water

---