

Yeast Protein Extraction Sample Buffer

Alan Cone

Abstract

Fresh Protein Extraction Sample Buffer.

Citation: Alan Cone Yeast Protein Extraction Sample Buffer. **protocols.io**

dx.doi.org/10.17504/protocols.io.ebmbak6

Published: 18 Dec 2015

Protocol

Step 1.

Add 2.4 mL 1 M Tris pH 6.8

 [AMOUNT](#)

2 ml Additional info:

Step 2.

Add 20 mL 50% Glycerol

 [AMOUNT](#)

20 ml Additional info:

 [REAGENTS](#)

 Glycerol [G5516](#) by Contributed by users

Step 3.

Add 8 mL 10% SDS

 [AMOUNT](#)

8 ml Additional info:

Step 4.

Add 4 mL of 1% BPB

 [AMOUNT](#)

4 ml Additional info:

Step 5.

Add 65.6 mL H₂O to make 100 mL of Buffer

 [AMOUNT](#)

66 ml Additional info: