

# Buffers for Chloroplast Isolation from Diatoms

P. Dreux Chappell, Bethany D. Jenkins

## Abstract

Guide for reagent/buffer preparation for protocol to separate chloroplasts from diatom cells using ammonium fluoride to permeate the silica frustule and a percoll gradient to separate the plastid from other cellular components.

**Citation:** P. Dreux Chappell, Bethany D. Jenkins Buffers for Chloroplast Isolation from Diatoms. **protocols.io**

dx.doi.org/10.17504/protocols.io.ea2bage

**Published:** 21 Apr 2016

## Guidelines

Buffers and reagents can be mixed up ahead of time and stored cold (refrigerated). Percoll gradients (layering of 40% and 80% percoll solutions) should be done the same day of the extraction shortly before use.

## Protocol

### Isolation Buffer (40 mL)

#### Step 1.

#### Isolation Buffer (40 mL)

1) Mix:

- 20 ml 1 M sorbitol
- 400 µl 0.6 M Na<sub>2</sub>-EDTA
- 200 µl 1 M MgCl<sub>2</sub>
- 400 µl 1 M KCl
- 40 µl 1 M MnCl<sub>2</sub>
- 2 ml 1 M HEPES-KOH pH = 8.0
- 16.96 ml Sterile H<sub>2</sub>O.

2) Add bovine serum albumen (BSA) 1% (w/v) just before use to subset of isolation buffer needed at the beginning of the procedure

### PEG-6000 Solution

#### Step 2.

#### PEG-6000 Solution

Mix:

- 6 g PEG

- 10 ml water

#### Percoll Solution (20 mL Stock)

##### Step 3.

##### Percoll Solution (20 mL Stock)

Mix:

- 19 ml Percoll
- 1 ml PEG – 6000 solution
- 0.2 g Ficoll
- 0.2 g BSA

#### Gradient Mixture (10 mL Stock)

##### Step 4.

##### Gradient Mixture (10 mL Stock)

Mix:

- 250  $\mu$ l 1 M HEPES-KOH pH = 8.0
- 1 ml 0.1 M EDTA
- 6.26 ml 1 M sorbitol
- 2.49 ml water

#### 3X Bottom Layer (80% Percoll)

##### Step 5.

##### 3X Bottom Layer (80% Percoll)

Mix:

- 5.7 ml Percoll Solution
- 1.01 ml Gradient Mixture

#### 3X Top Layer (40% Percoll)

##### Step 6.

##### 3X Top Layer (40% Percoll)

Mix:

- 5.04 ml Percoll Solution
- 6.96 ml Gradient Mixture

#### 10M Ammonium Fluoride (NH<sub>4</sub>F)

##### Step 7.

##### 10M Ammonium Fluoride (NH<sub>4</sub>F)

Mix:

- 3.7 g NH<sub>4</sub>F
- 10 ml H<sub>2</sub>O