



Sep 25, 2019

DNA Gel Electrophoresis

NUS iGEM¹¹National University of Singapore

1

Works for me

dx.doi.org/10.17504/protocols.io.7pshmne



NUS iGEM

National University of Singapore



STEPS MATERIALS

NAME ▾



CATALOG # ▾

VENDOR ▾

DNA Gel Loading Dye (6X)

R0611

Thermo Fisher Scientific



- 1 Load  10 µl of 6X loading dye into  50 µl of P product



DNA Gel Loading Dye (6X)

by Thermo Fisher Scientific

Catalog #: R0611

- 2 Mix and load entire volume into 1% agarose gel
- 3 Load  10 µl of 1kb or 100bp ladder
- 4 Run gel at 120V for  00:20:00
- 5 Visualize band under blue light transilluminator



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited