

# Quick staining procedure of nuclei in Euplotes using DAPI

Version 5

Rachele Cesaroni

## Abstract

**Citation:** Rachele Cesaroni Quick staining procedure of nuclei in Euplotes using DAPI. **protocols.io**

dx.doi.org/10.17504/protocols.io.hjfb4jn

**Published:** 01 Apr 2017

## Protocol

### Step 1.

Mix concentrated Euplotes cells together with Ethanol 70% in a ratio of 1:1.

### NOTES

**Rachele Cesaroni** 01 Apr 2017

Better to have a completely starved Euplotes culture to avoid autofluorescence from bacteria/algae.

### Step 2.

Add DAPI (0.01 mg/ml) to the mix in a ratio of 1:10, and stain for 15 minutes at room temperature.

### Step 3.

Observe cells by fluorescence microscopy.