

Co-Incubation protocol for transforming heterotrophic dinoflagellates (e.g. *Oxyrrhis marina*) Version 2

Brittany Sprecher, Lu Wang, Huan Zhang, and Senjie Lin

Abstract

Citation: Brittany Sprecher, Lu Wang, Huan Zhang, and Senjie Lin Co-Incubation protocol for transforming heterotrophic dinoflagellates (e.g. *Oxyrrhis marina*). **protocols.io**
dx.doi.org/10.17504/protocols.io.hmzb476

Published: 10 Apr 2017

Protocol

Prepare vector and transform into *E. coli*

Step 1.

Prepare a plasmid carrying proper promoter, selection marker, and reporter gene. Transform this plasmid into *coli* cells using standard heat shock protocol.

Grow transformed *E. coli*

Step 2.

Grow *coli* cells at a volume equal to 1/10 the desired volume of dinoflagellate transformation culture until the *E. coli* culture reaches an $OD_{600}=0.7$.

Harvest transformed *E. coli* cells

Step 3.

Centrifuge the cells at 3,000g for 5 mins to pellet the *E. coli*.

Mix transformed *E. coli* cells with prey alga

Step 4.

Discard the supernatant and resuspend the cells to twice the volume with either the favored prey of the species (*Dunaliella*) for transformation or in L1 media seawater.

Mix transformed *E. coli* cells and prey alga with dinoflagella

Step 5.

Add the *coli*-prey or -seawater mix to the dinoflagellate for transformation.

Co-incubation

Step 6.

After incubation at normal culture condition for six hours, add LB with a volume equal to 1/10 the total

volume. This step is not universal for all heterotrophic dinoflagellates so should be tested with every new species.

Add fresh medium and selection agent

Step 7.

After 24 hours add the same volume of new media and add the selection reagent to the media. Depending on the species it may be better to add the selection reagent after 42 hours.

Add fresh medium and selection agent

Step 8.

After 42 hours observe the culture under blue light.