

# western blot

Ning Chen

## Abstract

**Citation:** Ning Chen western blot. **protocols.io**

dx.doi.org/10.17504/protocols.io.iy5cfy6

**Published:** 18 Jul 2017

## Protocol

### Step 1.

Preparation of lysate from endometrial tissue—us T-PER Tissue Protein Extraction Reagent—Remove a small volume of lysate to perform a protein quantification assay. Determine the protein concentration for each tissue lysate—us BCA Protein Assay Kit.

### Step 2.

Load equal amounts of protein into the wells of the SDS -PAGE gel, along with molecular weight marker—60 ug total proteins per pore—Run the gel for 90min at 100 V—Gel percentage —separation gel 10%—spacer gel 5%.

### Step 3.

Activate PVDF with methanol for 1 min and rinse with transfer buffer before preparing the stack—Run at 100 V for 90min.

### Step 4.

Wash the membrane in three washes of TBST—5 percent skimmed milk powder—5 min each.

### Step 5.

Block the membrane for overnight at 4°C using blocking buffer. Wash the membrane in four washes of TBST, 5 min each.

### Step 6.

Incubate the membrane with the recommended dilution of conjugated secondary antibody in blocking buffer at room temperature for 1 h. Wash the membrane in five washes of TBST, 5 min each—

### Step 7.

Acquire image using normal image scanning methods for colorimetric detection.□

**Step 8.**

**Step 9.**