

# Protein secretion efficiency of construct in *Chlamydomonas reinhardtii* with fluorescent protein

João Vitor Molino

## Abstract

This protocol describe how to assess secretion efficiency of different constructs in *Chlamydomonas reinhardtii*, using a fluorescent protein. Colonies are picked from transformation plate and inoculated on TAP media for growth. A high number of colonies is picked to avoid bias on colony picking.

**Citation:** João Vitor Molino Protein secretion efficiency of construct in *Chlamydomonas reinhardtii* with fluorescent protein. **protocols.io**

[dx.doi.org/10.17504/protocols.io.kfnctme](https://dx.doi.org/10.17504/protocols.io.kfnctme)

**Published:** 01 Nov 2017

## Guidelines

Cell cultures in 96 well plates should ideally be cultivated at high orbital shaking (>800RPM) to avoid suboptimal growth.

## Before start

- Separate 96 well plates
- Prepare sufficient media
- Separate gas exchangeable seal

## Protocol

### Cell growth

#### Step 1.

1. Pipette 500  $\mu$ L of TAP media to each well of a 96 deep-well plate
2. Pick single colonies with autoclavated pipette tips or teeth cleaning twigs. Add to the wells
3. Cover plate with a gas exchangeable seal (Breathe-Easy®, Sigma-Aldrich®)
4. Place plates on a rotary shaker (150 RPM) under constant illumination (50  $\mu$ mol photons/m<sup>2</sup>s), at 25 °C for 7 days

#### AMOUNT

500  $\mu$ L Additional info: TAP media per well

## TEMPERATURE

25 °C Additional info: growth temperature

## NOTES

**João Vitor Molino** 01 Nov 2017

Low orbital shaker speed are sufficient for cell growth. Nevertheless, microtiter plate shakers can be used and enhance cell growth.

## Cell and supernatant separation

### Step 2.

1. Remove seal and transfer 100 µL of culture to fluorescence reading in a clear bottom 96-well plate (Corning Costar, Tewksbury, MA, USA).
2. Centrifuge the 96 deep-well plate at 3000xg for 10 min
3. Transfer 100 uL of supernatant to fluorescence reading in a clear bottom 96-well plate (Corning Costar, Tewksbury, MA, USA).

## AMOUNT

100 µl Additional info: sample size for each material

## Fluorescence reading

### Step 3.

1. Read fluorescence in the bottom of the plates at fluorescent protein specific wavelenghts in a plate reader
2. Ex: mCherry fluorescent protein
  - Excitation: 575/9 nm
  - Emission: 608/20 nm
  - Gain: 200
3. Read fluorescence in the bottom of the plates at chlorophyll specific wavelenghts in a plate reader
4. Ex: Chlorophyll
  - Excitation: 440/9 nm
  - Emission: 680/20 nm
  - Gain: 100

## NOTES

**João Vitor Molino** 24 Oct 2017

\*Compare results of the whole culture fluorescence and supernatant fluorescence to estimate secretion efficiency.

**\*\*Check chlorophyll fluorescence to surveil for contaminations or cell growth failure.**

## Absorbance reading

### Step 4.

1. Read cell density by absorbance at 750 nm

#### 📌 NOTES

**João Vitor Molino** 01 Nov 2017

Cell density can be used to normalize protein fluorescence results, and obtain a estimate of each colony productivity.