



Oct 15, 2019

Quantifying Biogenic Silica (bSi) Deposition Rates Adapted Method & Fluorescence Reading (PDMPO) via Fluorometer ↗

Brittany Zepernick¹, Matthew Saxton², Steven Wilhelm³

¹University of Tennessee, Knoxville, ²Miami University of Ohio, ³The University of Tennessee, Knoxville

1 Works for me dx.doi.org/10.17504/protocols.io.735hqq6

The Aquatic Microbial Ecology Research Group - AMERG (The Buchan, Zinser and Wilhelm labs)

Great Lakes Center for Fresh Waters and Human Health



Brittany N. Zepernick
University of Tennessee, Knoxville



ABSTRACT

This method can be used to assess and quantify the rate of silica deposition (bSi) over time in diatoms to determine their rate of frustule synthesis. This protocol has been adapted for the processing of both cultures as well as environmental samples when inoculated with PDMPO [2-(4-pyridyl)-5-((4-(2dimethylaminoethylaminocarbonyl)methoxy)phenyl)oxazole], which is a fluorescent dye that is co-deposited with silica during frustule biosynthesis in a 3230:1 Si:PDMPO (mol:mol) ratio.

EXTERNAL LINK

<http://wilhelmlab.utk.edu/>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Brittany N. Zepernick¹, Matthew A. Saxton², Steven W. Wilhelm¹ University of Tennessee Knoxville¹, Miami University² Adapted from: (Saxton et al, 2012), (LeBlanc and Hutchins 2005) Original Methods Leblanc, K., & Hutchins, D. A. (2005). New applications of a biogenic silica deposition fluorophore in the study of oceanic diatoms. *Limnology and Oceanography: Methods*, 3(10), 462-476. Saxton, M. A., D'souza, N. A., Bourbonniere, R. A., McKay, R. M. L., & Wilhelm, S. W. (2012). Seasonal Si: C ratios in Lake Erie diatoms—evidence of an active winter diatom community. *Journal of Great Lakes Research*, 38(2), 206-211.

MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

LysoSensor™ Yellow/Blue DND-160 - Special Packaging

L7545

Fisher Scientific

MATERIALS TEXT

Utilize a TD-700 Laboratory Fluorometer linked to a computer monitor with the downloaded Turner software to obtain fluorescence readings of the samples

SAFETY WARNINGS

See SDS (Safety Data Sheet) for hazards and safety warnings.

Collection protocols



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited