



Murashige and Skoog (MS) agar

Forked from [Murashige and Skoog \(MS\) medium](#)

[Steven Burgess](#)¹

¹University of Illinois at Urbana-Champaign

[dx.doi.org/10.17504/protocols.io.xd8fi9w](https://doi.org/10.17504/protocols.io.xd8fi9w)

Plantae



Steven Burgess

University of Illinois at Urbana-Champaign



ABSTRACT

Murashige and Skoog medium (or **MSO** or **MS0** (*MS-zero*)) is a [plant growth medium](#) used in the laboratories for cultivation of plant [cell culture](#). MSO was invented by plant scientists [Toshio Murashige](#) and [Folke K. Skoog](#) in 1962 during Murashige's search for a new [plant growth regulator](#). A number behind the letters MS is used to indicate the sucrose concentration of the medium. For example, MS0 contains no sucrose and MS20 contains 20 g/l sucrose. Along with its modifications, it is the most commonly used medium in [plant tissue culture](#) experiments in laboratory.^[1]

Source:

http://cshprotocols.cshlp.org/content/2010/2/pdb.rec12142.full?text_only=true

EXTERNAL LINK

https://en.wikipedia.org/wiki/Murashige_and_Skoog_medium

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

Derivatives of MS medium can be made through the addition of sucrose, both 1% (w/v) and 2% (w/v) are commonly used.

- 1 Add Murashige and Skoog Basal Salt medium to a 1L flask **4.33 g**
- 2 Add dH₂O **800 ml**
- 3 Adjust the pH to 5.7 using 2 N Potassium hydroxide KOH
- 4 Add dH₂O up to 1L
- 5 Add **7 g** of Bacto Agar to flask and autoclave



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

