



Plasma-derived exosomal preparation and isolation of RNA

PLOS One

Fernanda Svedman¹, Warangkana Lohcharoenkal¹, Matteo Bottai¹, Suzanne Egyhazi Brage¹, Enikö Sonkoly¹, Johan Hansson¹, Andor Pivarcsi¹, Hanna Eriksson¹

¹Karolinska Institute Stockholm

dx.doi.org/10.17504/protocols.io.tuienue



Fernanda Svedman



EXTERNAL LINK

https://doi.org/10.1371/journal.pone.0206942

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Svedman FC, Lohcharoenkal W, Bottai M, Brage SE, Sonkoly E, Hansson J, Pivarcsi A, Eriksson H (2018) Extracellular microvesicle microRNAs as predictive biomarkers for targeted therapy in metastastic cutaneous malignant melanoma. PLoS ONE 13(11): e0206942. doi: 10.1371/journal.pone.0206942

PROTOCOL STATUS

Working

Preparation of plasma-derived EVs and isolation of RNA

① Section A._3.png

@ Section A._2.png

MiRCURY Exosome Isolation Kit-Serum and Plasma (Exiqon, Denmark) was used for isolation of plasmatic EVs. Briefly, 600 ml of plasma was treated with thrombin for 5 minutes at room temperature to remove clotting factors. The supernatant was mixed with 200 ml of precipitation buffer and incubated at 4°C for overnight and centrifuged (10,000 g for 30 min at room temperature). The pellet was vortexed for 5-15 minutes in 270 ml resuspension buffer until completely re-suspended. The EV miRNA was extracted using miRCURY RNA isolation kit-biofluids (Exiqon, Denmark) according to the manufacturerÂ's instruction. RNA spike-in template mixture (Exiqon, Denmark) was added as a quality control of the downstream PCR analysis.

(I) Section A. Standard protocol for RNA isolation_1.png

(a) Section D. Exosome isolation.png

2

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited