



May 17,
2019

Working

Yale - Blood Glucose [↗](#)

John Stack¹, Gary Cline¹

¹Yale University

[dx.doi.org/10.17504/protocols.io.y2kfyw](https://doi.org/10.17504/protocols.io.y2kfyw)

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

Summary:

Procedure used to measure glucose concentrations. Glucose is measured by the enzymatically coupled reactions of hexokinase and glucose-6-P dehydrogenase. The rate of NADH formation is monitored by the change in absorbance at 340 nm.

EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=202&docType=Protocol>

MATERIALS

NAME	CATALOG #	VENDOR
Glucose Reagent 1	R84682	Prolabs(cliniqa)
Glucose Reagent 2	R84682	Prolabs(cliniqa)
Multi Analyte Calibrator	R60010	Prolabs(cliniqa)
Assayed Control Serum 1	R83082	Prolabs(cliniqa)
Assayed Control Serum 2	R83083	Prolabs(cliniqa)

MATERIALS TEXT

Reagent Preparation:

Glucose Reagent 1: As supplied by Vendor

Glucose Reagent 2: As supplied by Vendor

Multi Analyte Calibrator: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Assayed Control Serum 1: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Assayed Control Serum 2: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Note: Author switched vendors and are now using a [Sekisui 235-60 reagent](#) for glucose assay. (There is now only 1 reagent instead of 2).

- 1 Calibrate Cobas for Glucose analysis by running a multi analyte standard and two control serum.

- 2 Sample Handling as performed by the Cobas Mira Plus.
 - a) Pipette 3 μ L of sample into cuvette.
 - b) Absorbance is measured at 340 nm.
 - c) Add 100 μ L of Glucose liquid reagent.
 - d) Mixture is incubated at 37°C for 10 minutes.
 - e) Absorbance is measured at 340 nm. Change in absorbance is calculated.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited