

Seydoux lab Cas9 preparation test: Tag *gtbp-1* with eGFP or mCherry using *dpy-10* co-CRISPR

Alexandre Paix, Andrew Folkmann, Dominique Rasoloson, Geraldine Seydoux

Abstract

This is a positive control experiment (to test protocol in your hands and/or activity of your home-made Cas9): tag *gtbp-1* with eGFP or mCherry using *dpy-10* co-CRISPR.. The protocol is from:

Paix A, Folkmann A, Rasoloson D, and Seydoux G (2015) [High Efficiency, Homology-Directed Genome Editing in *Caenorhabditis elegans* Using CRISPR/Cas9 Ribonucleoprotein Complexes](#). Genetics genetics.115.179382; Early online July 17, 2015. doi:10.1534/genetics.115.179382

Please see the [full manuscript](#) for detail.

Citation: Alexandre Paix, Andrew Folkmann, Dominique Rasoloson, Geraldine Seydoux Seydoux lab Cas9 preparation test: Tag *gtbp-1* with eGFP or mCherry using *dpy-10* co-CRISPR. **protocols.io**
dx.doi.org/10.17504/protocols.io.duj6um

Published: 06 Oct 2015

Guidelines

Main protocol: [Direct delivery CRISPR-HDR editing protocol for *C. elegans*](#)

Reagents for *gtbp-1* eGFP and mCherry tagging

crRNA *gtbp-1* Ct (#728 in Table S3): CCACGAGGUGGUAUGCGCAG + GUUUUAGAGCUAUGCUGUUUUG

Primers pairs to generate template (from eGFP containing-plasmid pAP682-1):

Forward primer (5' to 3', lower case indicating the homology arm sequence) for eGFP insertion:
ggttcgggtggtgctccacgagtggtatgctgGTGAGTAAAGGAGAAGAAC

Reverse primer (5' to 3', lower case indicating the homology arm sequence) for eGFP insertion:
cttctaattttgtcccgcatTTTggaaaccgctCTTGACAGCTCGTCCATGCC

Primers pairs to generate template (from mCherry containing-plasmid pAP582-1):

Forward primer (5' to 3', lower case indicating the homology arm sequence) for mCherry insertion:
ggttcgggtggtgctccacgagtggtatgctgGTCTCAAAGGGTGAAGAAGATAAC

Reverse primer (5' to 3' lower case indicating the homology arm sequence) for mCherry insertion:
cttctaattttgtcccgcatTTTggaaaccgctCTTATACAATTCATCCATGCC

PCR condition: use annealing temperature of 63 °C, elongation step of 45s, pAP682-1 (eGFP) or

pAP582-1

(mCherry) as a plasmid template. Do 8 PCR reactions and pool them (400µl total), purify on minelute column (see [Repair template synthesis](#)).

Protocol

Step 1.

Prepare the repair PCR template as indicated in **Reagents for gtbp-1 eGFP and mCherry tagging** (see guidelines) and [Repair template synthesis](#).

Step 2.

Make the injection mix as indicated in [injection mixes, part I](#). Use the crRNA gtbp-1 Ct (#728 in Table S3).

Step 3.

Inject 15-20 young adult N2 worms and recover as described in [Worm recovery and handling](#).

Step 4.

Pool the recovered worms (P0s) on one plate and incubate for 22-23h at 20 °C (day 1).

 **DURATION**

22:00:00

Step 5.

Clone the P0s to individual OP50 plates and incubate at 20 °C for 4-5 days (day 2).

Step 6.

When the F1s reach the adult stage, check for Rollers.

 **NOTES**

Alexandre Paix 18 Sep 2015

At least 3 P0s should give a high number of Rollers (>15)

Alexandre Paix 18 Sep 2015

Note that we do not count/examine Dumpy F1s since these are homozygous edits at the *dpy-10* locus.

Step 7.

From those “jackpot broods”, screen the Rollers for fluorescent protein expression as described in [Screening for edit of interest](#).

 **NOTES**

Alexandre Paix 18 Sep 2015

50% or more of the Rollers (at day 2) should be positive for fluorescence.