

Revisited *Thalassiosira pseudonana* (Tp) conjugation protocol enables fusion protein delivery to Tp frustule version 3

Jernej Turnsek

Abstract

This revisited *Thalassiosira pseudonana* conjugation protocol builds on the following two protocols: [Conjugation of *Thalassiosira pseudonana*](https://doi.org/10.17504/protocols.io.f55bq86) (DOI: [dx.doi.org/10.17504/protocols.io.f55bq86](https://doi.org/10.17504/protocols.io.f55bq86)) & [Pour plating of *Thalassiosira pseudonana* \(Tp\)](https://doi.org/10.17504/protocols.io.jfncjme) (DOI: [dx.doi.org/10.17504/protocols.io.jfncjme](https://doi.org/10.17504/protocols.io.jfncjme)). Readers are strongly encouraged to read both before using this one. Gene-EGFP fusion depicted in the 'Expected Results' section was expressed under constitutive Tp-derived fcp promoter and terminator using pTpPuc3 episome backbone described in Karas *et al.* (2015).

Citation: Jernej Turnsek Revisited *Thalassiosira pseudonana* (Tp) conjugation protocol enables fusion protein delivery to Tp frustule. **protocols.io**

[dx.doi.org/10.17504/protocols.io.nbzdap6](https://doi.org/10.17504/protocols.io.nbzdap6)

Published: 18 Feb 2018

Protocol

Growth and preparation of *E. coli* (donor) cells

Step 1.

Pick bacterial colonies from Gent20 (pTA-Mob) + Kan50 (pTpPuc3 episome) plates and inoculate 3 mL LB medium. You can use pre-made glycerol stocks instead. Shake overnight at 37 °C.

Growth and preparation of *E. coli* (donor) cells

Step 2.

Passage 1:200 in 50 mL LB medium. Shake at 37 °C until OD₆₀₀ reaches 0.3–0.4.

Growth and preparation of *E. coli* (donor) cells

Step 3.

Centrifuge cells 10 min / 4000 rpm / 10 °C.

Growth and preparation of *E. coli* (donor) cells

Step 4.

Decant all supernatant and resuspend pellet in 267 µL SOC medium.

NOTES

Jernej Turnsek 18 Feb 2018

After concentrating *E. coli* cells, it is important to remove all excess supernatant as it contains antibiotics that could potentially alter the conjugation reaction.

Growth and preparation of *T. pseudonana* (recipient) cells

Step 5.

Determine cell density of your *Tp* culture. 1.3×10^7 cells/mL "stock" culture was used in this protocol.

📌 NOTES

Jernej Turnsek 17 Feb 2018

'Tp L1' was buffered with 10 mM Tris-HCl and pH was adjusted to 8. This allows longer maintenance of *Tp* in culture without crashing. It also induces sexual reproduction, possibly due to excess nitrate in 'Tp L1' and extra nitrogen from Tris. Auxospores can be readily observed in this buffered 'Tp L1'. Readers are referred to [this study](#) for further information.

'Tp L1' preparation details are available in [this protocol](#).

Growth and preparation of *T. pseudonana* (recipient) cells

Step 6.

Centrifuge 2×10^8 cells 10 min / 4000 rpm / 10 °C. Add a few mL to your calculation due to cell death and losses during centrifugation (e.g. in this protocols 20 instead of calculated 15.4 mL were spun down).

Growth and preparation of *T. pseudonana* (recipient) cells

Step 7.

Resuspend cells in 1 mL "Tp L1".

Conjugation

Step 8.

Add 200 μ L (4×10^7 cells) *Tp* cell suspension to 1.5 mL eptube.

Conjugation

Step 9.

Add 200 μ L *E. coli* cell suspension. Pipette up and down a few (5) times.

Conjugation

Step 10.

Plate the mix (i.e. all 400 μ L) on 1/2 'Tp L1' 5% LB 1% agar plate.

📌 NOTES

Jernej Turnsek 17 Feb 2018

Plates have to be pre-dried. Leave them half lid open in the hood for 90 min before plating.

Conjugation

Step 11.

Incubate plates 90 min at 30 °C in dark (i.e. where you normally grow yeast).

Conjugation

Step 12.

Incubate plates additional 4 hrs at 18 °C under constant light (i.e. where you normally grow *Tp*).

Conjugation

Step 13.

Scrape co-culture with 1 mL 'Tp L1' and transfer to 1.5 mL eptube. Use [L spreader](#). Expect 600–700 µL material.

Selection

Step 14.

Plate 50 and 500 µL using [the pour plating method](#) (**crucial**). Use 100 µg/mL nourseothricin sulfate for selection in pour plates.

📌 NOTES

Jernej Turnsek 17 Feb 2018

1. Regular non-superclean agar may work, but this hasn't been tested (yet).
2. Some WT *Tp* colonies were observed in/on control plates which could be due to antibiotic depletion (2×10^7 cells were plated). Control plates contained 100 µg/mL nourseothricin sulfate.

Brief summary of results

Step 15.

Dozens of colonies emerged in pour plates after 2 weeks and 300 µL 'minicultures' were started from single colonies (200 µg/mL, not 100, nourseothricin sulfate to additionally increase selection pressure). Gene-specific primers were used for genotyping and Leica TCS SP5 for EGFP signal screening. Scale bars below: 25 µm (top row), 5 µm. Expected cell wall localization pattern was observed in select lines.

EXPECTED RESULTS

