



Mar 22,  
2019

Working

## Removal of genomic DNA from RNA preparations (Thermo Scientific )- (M4455 Version)

Forked from [Removal of genomic DNA from RNA preparations \(Thermo Scientific \)](#)

Maximilian Dietsch<sup>1</sup>

<sup>1</sup>Heinrich-Heine Universität Düsseldorf

[dx.doi.org/10.17504/protocols.io.ze2f3ge](https://doi.org/10.17504/protocols.io.ze2f3ge)

M4455 - Synthetische Biologie und Biotechnologie



Sebastian Triesch

Institute for Synthetic Microbiology, HHU Düsseldorf



### ABSTRACT

Removal of genomic DNA from RNA preparations

### PROTOCOL STATUS

#### Working

We use this protocol in our group and it is working

### GUIDELINES

always waer gloves and work on ice

### SAFETY WARNINGS

1	<b>Add to an RNase free tube:</b>	
	RNA	1 µg
	10X reaction buffer with MgCl <sub>2</sub>	1 µl
	DNase I, RNase-free	1 µl (1U)
	Water	to 10 µl

2 Incubate at 37 °C for 30 min

37 °C

00:30:00

3 Add EDTA, Water and PCI and vortex thoroughly.

1 µl EDTA

80 µl Water

100 µl PCI (phenol chloroform isoamyl alcohol)

4 Centrifuge for 10 min at 10000 rpm and 4 °C

5 transfer the upper phase into a fresh tube and add 3 volumen EtOH/ 3M Natrumacetat (30:1, ph 5.2 )

6 precipitate RNA over night at -20 °C

7 Centrifuge 30 min at 13000 rpm and 4 ° C

 4 °C

 00:30:00

8 Discard supernatant and wash pellet with 75% EtOH ( do not resuspend the pellet)

9 Centrifuge 10 min at 13000 rpm and 4 °C

 4 °C

 00:10:00

10 [go to step #8 repeat washing step](#)

[go to step #9 Centrifuge](#)

11 Discard supernatant and dry pellet for 10 - 15 min

12 resuspend pellet with 30 µl H<sub>2</sub>O



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited