

# Freezing and unfreezing unicellular protists using Mr Frosty container

Maria Rubio-Brotons

## Abstract

**Citation:** Maria Rubio-Brotons Freezing and unfreezing unicellular protists using Mr Frosty container. **protocols.io**  
dx.doi.org/10.17504/protocols.io.ggybtxw

**Published:** 22 Nov 2016

## Protocol

Collect your culture in 15 mL falcon tube and centrifuge at 1000 g for 5 minutes

### Step 1.

Discard medium and resuspend in 10 % DMSO growing medium

### Step 2.

Use growing medium according to your organism, in our case is marine broth.

We usually make 10 mL final volume (1mL DMSO into 9mL growing medium).

Aliquot 1 mL of your culture into criogenic vials

### Step 3.

## Nalgene<sup>®</sup> cryogenic vials

Place your criovials into Mr Frosty containing 100% isopropanol

### Step 4.

Nalgene<sup>®</sup> Mr. Frosty<sup>®</sup> Cryo 1°C Freezing Containers

Store Mr Frosty container at -80°C for at least 4 h

### Step 5.

We usually keep it for 24 hours

Move your criovials into a long-term storage container at -80°C

### Step 6.

Unfreezing your culture

### Step 7.

Take your criovials out of -80°C storage and place them into 35°C water bath for 2-3 minutes

**Step 8.**

When you see they start thawing take them out of the water bath, usually is after 2-3 minutes

Transfer into fresh growing medium

**Step 9.**

We usually have 1mL frozen culture, transfer this 1 mL into 7 mL fresh growing medium flask

Monitor survival rates within the next 24h

**Step 10.**

After 2-3 days unfreezing, pellet your cells and resuspend with fresh growing medium

**Step 11.**

Centrifuge 1000 g for 5 minutes