

DNA/RNA Extractions from Sterivex Filters

Matthew Sullivan

Abstract

Citation: Matthew Sullivan DNA/RNA Extractions from Sterivex Filters. **protocols.io**

dx.doi.org/10.17504/protocols.io.c4bysm

Published: 21 Jan 2016

Guidelines

2.0 mL Lysis Buffer is added to Sterivex immediately after filtration. Filters are stored at -20°C or -80°C.

Lysis Buffer is:

Final Concentration For 20 ml:

40 mM EDTA 1.6 ml of 0.5 M EDTA

50 mM Tris (pH 8.3) 1.0 of 1 M Tris (pH 8.3)

0.73 M Sucrose 5.13g sucrose

* Note: Must use Phenol:Chloroform:IAA (25:24:1) pH 8.0! Add buffer supplied, shake well, allow phases to separate and check pH of buffer using pH paper. Add hydroxyquinoline to a final concentration of 0.1% (100 mg hydroxyquinoline per 100ml phenol:chloroform:IAA).

Before start

Prepare lysis buffer as described in guidelines.

Protocol

Step 1.

Thaw sterivex filters on ice.

Step 2.

Add 40µL lysozyme solution (add 2 mg lysozyme to 40 µl Lysis Buffer) to filter.

 **AMOUNT**

40 µl Additional info:

Step 3.

Incubate at 37°C while rotating for 45 minutes.

 **DURATION**

00:45:00

Step 4.

Add 100µL Proteinase K solution (1mg Proteinase K in 100 µl Lysis Buffer) and 100µL 20% SDS.

AMOUNT

100 µl Additional info:

Step 5.

Incubate at 55°C while rotating for 2 hours.

DURATION

02:00:00

Step 6.

Transfer lysate to a 15 mL conical tube using a sterile 3 mL syringe.

Step 7.

Add 1 mL Lysis Buffer to filter and wash at 55°C for 15 minutes. Pool with above lysate.

AMOUNT

1 ml Additional info:

DURATION

00:15:00

PROTOCOL

. [Lysis Buffer \(20 mL\)](#)

CONTACT: [Celina Gomez](#)

Step 7.1.

1.6 mL of 0.5 M EDTA

REAGENTS

EDTA disodium dihydrate [AB1011793](#) by [Abblis](#)

Step 7.2.

1.0 mL of 1 M Tris (pH 8.3)

REAGENTS

 Tris [RP-T60040](#) by [P212121](#)

Step 7.3.

5.13 g Sucrose

REAGENTS

 Sucrose [View](#) by [P212121](#)

Step 8.

Add 3 mL Phenol:Chloroform:IAA (25:24:1; pH 8.0).

AMOUNT

3 ml Additional info:

NOTES

VERVE Team 18 Jun 2015

Must use Phenol:Chloroform:IAA (25:24:1) at pH 8.0. Add buffer supplied, shake well, allow the phases to separate and check pH of buffer using pH paper. Add hydroxyquinoline to a final concentration of 0.1% (100 mg hydroxyquinoline per 100 mL Phenol:Chloroform:IAA).

Step 9.

Vortex for 10 seconds.

 DURATION

00:00:10

Step 10.

Spin for 5 minutes at 2500g (speed 8).

 DURATION

00:05:00

Step 11.

Carefully transfer aqueous phase to a new 15 mL conical tube.

Step 12.

Add 3 mL Chloroform:IAA (24:1)

 AMOUNT

3 ml Additional info:

 REAGENTS

Chloroform:IAA [C0549-1PT](#) by [Abblis](#)

Step 13.

Vortex for 10 seconds.

 DURATION

00:00:10

Step 14.

Spin for 5 minutes at 2500g (speed 8).

 DURATION

00:05:00

Step 15.

Carefully transfer aqueous phase to a centricon 100.

Step 16.

All of the remaining volume should be transferred to an epi tube.

Step 17.

Spin Centricon at 1000xg (speed4) for 20 minutes.

 DURATION

00:20:00

Step 18.

Add remaining volume from epi tube to centricon.

Step 19.

Spin until only 100µL - 500µL of aqueous phase is left in Centricon.

Step 20.

Add 1 mL TE Buffer

 AMOUNT

1 ml Additional info:

Step 21.

Repeat step 20.

Step 22.

Spin until only 100 μ L - 150 μ L left in Centricon.

Step 23.

Carefully remove liquid without damaging the membrane.

Step 24.

Wash membrane well with 40 μ L of TE Buffer

 AMOUNT

40 μ L Additional info:

Step 25.

Pool the membrane in the same epi tube. Note the total volume in the epi tube.