



Jun 11,
2019

Prep for FACS on blood samples

Grace Burgin¹, Inbal Benhar¹

¹Broad Institute

Working

[dx.doi.org/10.17504/protocols.io.32dgqa6](https://doi.org/10.17504/protocols.io.32dgqa6)



Grace Burgin



ABSTRACT

To prepare blood samples for flow

- 1 Collect blood into EDTA coated tubes
- 2 Transfer 50ul into clean tube - include 1 unstained
***mix from several mice to have 50ul for unstained ctrl tube
- 3 Add 1ml ACK, incubate at RT for 7 min, centrifuge 5' 350g
- 4 Resuspend pellet in FACS buffer + antibody mix
- 5 Incubate in dark for 15' on ice
- 6 Wash (centrifuge) in FACS buffer, 5' 350g
- 7 Resuspend in 300ul FACS Buffer + calcein blue



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited