

# Recovery of DNA from Low Melting Point Agarose

Irina Agarkova

## Abstract

**Citation:** Irina Agarkova Recovery of DNA from Low Melting Point Agarose. **protocols.io**

[dx.doi.org/10.17504/protocols.io.et2beqe](https://dx.doi.org/10.17504/protocols.io.et2beqe)

**Published:** 13 Jun 2016

## Protocol

### Step 1.

Low melting point agarose gels should be poured and allowed to set up in the cold room, but electrophoresed at room temperature.

### Step 2.

The gels should be stained at 4°C and kept cold until ready to cut out the DNA bands from the gel.

### 📌 NOTES

**Irina Agarkova** 05 Apr 2016

The gels will be fragile, so handle gently.

### Step 3.

Cut out the DNA bands from the gel with new razor blades. Use a new razor blade for each band.

### Step 4.

Place the agarose pieces into sterile capped tubes (the size of the tubes depends on the volume of the agarose involved).

### Step 5.

Add 1X TE buffer to the tubes of agarose (as small a volume as possible, usually about 5X the volume of the gel pieces).

### Step 6.

Heat the tubes at 65°C for 10-15 min.

### 🕒 DURATION

00:15:00

### Step 7.

Transfer the tubes to 37°C.

### Step 8.

Add an equal volume of buffer-saturated phenol. The phenol should have been warmed to 37°C.

### Step 9.

Mix well.

### Step 10.

Centrifuge the tubes at room temperature to separate the phases. Centrifuge at 10,000 rpm, 10 min in the Sorvall or 5 min in the microfuge.

### 🕒 DURATION

00:10:00

**Step 11.**

Transfer the upper aqueous layers to clean tubes.

**Step 12.**

Add an equal volume of phenol:CHCl<sub>3</sub>:Isoamyl alcohol (25:24:1) to the tubes.

**Step 13.**

Mix well and centrifuge at 10,000 rpm, 10 min in the Sorvall or 5 min in the microfuge.

 DURATION

00:10:00

**Step 14.**

Re-extract the aqueous layer 2X with CHCl<sub>3</sub>:Isoamyl alcohol (24:1) in the centrifuge as before.

**Step 15.**

Add 3 M Na acetate to a final concentration of 0.3 M and precipitate the DNAs with 2X volumes of 100% EtOH at -20°C overnight.

 DURATION

18:00:00

**Step 16.**

Centrifuge the tubes to pellet the DNAs. Centrifuge in the Sorvall at 10,000 rpm, 10 min, 4°C or 10-15 min at 4°C in the microfuge

 DURATION

00:10:00

**Step 17.**

Wash the DNAs 1X with 70% EtOH in the centrifuge as before and dry the DNA pellets in the vacuum desiccator briefly (10-15 min) to remove the EtOH.

 DURATION

00:15:00

**Step 18.**

Resuspend the DNAs with a small volume of 1X TE buffer.