

Protein extraction from Aiptasia

OSU Weis Lab

Abstract

This protocol, developed in 2011 by Angela Poole while in the Weis lab, is a quick prep for protein extraction from Aiptasia adults.

Citation: OSU Weis Lab Protein extraction from Aiptasia. protocols.io

dx.doi.org/10.17504/protocols.io.rifd4bn

Published: 23 Jul 2018

Before start

RIPA extraction buffer:

- 100mM Tris, pH 7.4
- 100mM NaCl
- 10mM EDTA, pH 8.0

*Bring to desired volume with Nanopure water and just before use, dissolve one cOmplete™, mini, protease inhibitor tablet in 10mL of RIPA Buffer. Be sure to let the tablet dissolve on its own- do not vortex!!

Materials



cOmplete™, Mini Protease Inhibitor Cocktail
11836153001 by Roche

Protocol

Step 1.

Prepare anemones by placing in the incubator in artificial seawater (ASW) 3-4 days prior to extraction. Be sure to change out the (ASW) each day to remove any residual debris.

Step 2.

Transfer 4-5 large anemones to a small tissue grinder on ice with 0.5-1 mL of chilled extraction buffer. Upon removal from their original container, anemones should be blotted on Kimwipes/weigh paper to remove as much water as possible.

Step 3.

Transfer homogenate to a 1.5 mL tube and centrifuge at 14,000xg for 15 minutes at 4°C. This step will pellet the dinoflagellates and cell debris, while the protein will be in the supernatant.

Step 4.

Remove the supernatant (be careful not to disturb the white lipid layer!) and place in a new tube.

Step 5.

Determine the concentration of your protein (we used the Bradford assay).

Step 6.

Aliquot your protein into 100 µl volumes and store in the -80°C freezer until further use. Freeze-thaw cycles should be avoided and protein can be kept for a short period in the fridge on ice.
