

# Transformation Protocol

## New England Biolabs

### Abstract

Quick Ligation products may be transformed by many different methods. The following protocol is recommended by New England Biolabs.

**Citation:** New England Biolabs Transformation Protocol. **protocols.io**

dx.doi.org/10.17504/protocols.io.cpxvpm

**Published:** 25 Jan 2015

## Protocol

### Step 1.

Thaw competent cells on ice

### Step 2.

Chill approximately 5 ng (2 µl) of the ligation mixture in a 1.5 ml microcentrifuge tube.

 **AMOUNT**

2 µl Additional info:

### Step 3.

Add 50 µl of competent cells to the DNA.

 **AMOUNT**

50 µl Additional info:

### Step 4.

Mix gently by pipetting up and down or flicking the tube 4–5 times to mix the cells and DNA. Do not vortex.

### Step 5.

Place the mixture on ice for 30 minutes. Do not mix.

 **DURATION**

00:30:00

### Step 6.

Heat shock at 42°C for 30 seconds. Do not mix.

 **DURATION**

00:00:30

 **ANNOTATIONS**

**New England Biolabs** 27 Jan 2015

For the duration and temperature of the heat shock step, please follow the recommendations provided by the competent cells' manufacturer.

### Step 7.

Add 950 µl of room temperature media to the tube.

 **AMOUNT**

950 µl Additional info:

## ■ ANNOTATIONS

**New England Biolabs** 27 Jan 2015

For the media to be used during the recovery period, please follow the recommendations provided by the competent cells' manufacturer.

### **Step 8.**

Place tube at 37°C for 60 minutes. Shake vigorously (250 rpm) or rotate.

🕒 **DURATION**

01:00:00

### **Step 9.**

Warm selection plates to 37°C.

### **Step 10.**

Spread 50–100 µl of the cells and ligation mixture onto the plates.

### **Step 11.**

Incubate overnight at 37°C

🕒 **DURATION**

15:00:00