

# Buffered glutaraldehyde

Li Deng

## Abstract

Formulation for a phosphate-buffered 20% glutaraldehyde stock solution.

**Citation:** Li Deng Buffered glutaraldehyde. **protocols.io**

dx.doi.org/10.17504/protocols.io.c7gzjv

**Published:** 30 Nov 2015

## Guidelines

Solution A

Solution B

$\text{Na}_2\text{HPO}_4$  1.414g in 50ml 0.02 $\mu\text{m}$  filtered mQ  $\text{NaH}_2\text{PO}_4$  1.058g in 50ml 0.02 $\mu\text{m}$  filtered mQ

## Protocol

### Step 1.

Prepare solution A and solution B as found in guidelines

### Step 2.

Prepare the phosphate buffere by adding 15.25 ml solution A and 9.75 ml solution B in 25 ml 0.02  $\mu\text{m}$  filtered mQ

#### 🔗 NOTES

**Bonnie Poulos** 24 Jun 2015

= 50 ml buffer

### Step 3.

Add 40 ml of 25% glutaraldehyde to 10 ml phosphate buffer for a 20% buffered glutaraldehyde solution

#### 🔗 NOTES

**Bonnie Poulos** 30 Nov 2015

This is most often used for transmission electron microscopy so purchase a TEM-grade glutaraldehyde (from a vendor such as Ted Pella) for preparing this buffered reagent

### Step 4.

Adjust pH to 7 with NaOH

### Step 5.

Store this 20% buffered glutaraldehyde at 4 °C in the dark