

Supernatant Transfer Assay and Transwell® Assay

Sandra L. Ross, Marika Sherman, Patricia L. McElroy, Julie A. Lofgren, Gordon Moody, Patrick A. Baeuerle, Angela Coxon, Tara Arvedson

Abstract

This protocol describes two methods for assessing whether soluble factors alone or a combination of soluble factors and T cells are responsible for bystander killing: (1) Supernatant transfer assays and (2) Transwell® assays.

Citation: Sandra L. Ross, Marika Sherman, Patricia L. McElroy, Julie A. Lofgren, Gordon Moody, Patrick A. Baeuerle, Angela Coxon, Tara Arvedson Supernatant Transfer Assay and Transwell® Assay. **protocols.io** dx.doi.org/10.17504/protocols.io.h3vb8n6

Published: 17 May 2017

Protocol

Overview

Step 1.

This protocol describes methods to assess the contribution of soluble factors and BiTE®-activated T cells in bystander killing. Two methods are used:

(1) Supernatant transfer assays from BiTE®- + T cell-treated target positive cells to target-negative cells and (2) Transwell® assays in which a membrane separates BiTE®- + T cell-treated target-positive cultures from target-negative cultures, restricting movement of T cells but not soluble factors.

These assays are modifications of the standard TDCC protocol. Please refer to TDCC protocol for cell treatments and assay set-up.

Materials

Step 2.

See also Materials list in TDCC protocol.

HTS Transwell-96 System, PET membrane Corning #3392 or #3380 1 µm

HTS Transwell-96 Well Plate, PC membrane Corning #3387 or #3388 5 µm

HTS Transwell-96 receiver plate, tissue culture treated Corning #3382 clear plate

HTS Transwell-96 receiver plate, tissue culture treated Corning #3583 black plate, including bottom

Supernatant transfer assay

Step 3.

For 96-well assays, each well contains a final volume of 200 µl. Plan for setting up cells in separate ViewPlates to allow for:

- Supernatant transferred directly (medium + cells)
- Supernatant clarified by centrifugation (medium only)
- No transfer control

1. Prepare target-positive cells:

- Detach adherent target cells, count and dilute cells (e.g. NUGC-4) in a mixture with T cells at desired E:T ratio (10:1 typically used) in 160 µl, e.g. 10:1 = 10,000 target cells and 100,000 T cells per well
- Add 40 µl per well of BiTE® dilution or assay medium to plates
- Add 160 µl per well target cells plus T cell mixture to plates
- Allow plates to equilibrate for 30 min. at room temperature
- Incubate at 37°C, 5% CO₂ for 48 hours

3. Prepare target-negative cell plates (e.g. SW620) on day of supernatant transfer:

- Detach adherent target-negative cells, count and dilute cells in a mixture with T cells at desired E:T ratio in 160 µl, e.g. 10:1 = 10,000 target cells and 100,000 T cells per well
- Add 160 µl target-negative cells plus T cell mixture to a fresh 96-well ViewPlate
- Allow plates to equilibrate for 30 min. at room temperature
- Incubate at 37°C, 5% CO₂ for 5 hours to allow cells to adhere.

4. Set up supernatant transfer assay:

- Remove media from target-negative cells by inverting plate and blot edges on paper towel
- Plate 1: Transfer directly 0.160 ml crude supernatant from target-positive NUGC-4 plate to plate containing target-negative SW620 cells
- Plate 2: Transfer carefully 0.160 ml clarified supernatant from target-positive NUGC-4 plate to plate containing target-negative SW620 cells. (NUGC-4 plate first centrifuged at 400 x g for 5 min.)
- Plate 3: No transfer control-add fresh BiTE® dilutions and T cells
- Incubate at 37°C, 5% CO₂ for 48 hours

5. Perform cytotoxicity assay (imaging or CellTiter-Glo® - see TDCC protocol).

Transwell® assay

Step 4.

For 96-well transwell assays, each well contains a final volume of 310 µl

- Bottom receiver plate: 200 µl target-negative (or target-positive for positive control) cells + 35 µl assay medium
- Top chamber insert: 35 µl 10X BiTE® dilutions + 40 µl target-positive cells/T cell mixture

(typically 10:1 E:T ratio)

1. Prepare target cells:

- Detach adherent target-positive cells, count and dilute to 10,000 cells/well, 200 µl/well. Also, dilute target cells in a mixture with T cells at desired E:T ratio in 40 µl, e.g. 10:1 = 10,000 target cells and 100,000 T cells per well

2. Set up assay:

- Add 35 µl assay medium to wells in bottom of receiver plate
- Add 200 µl target-negative (or target-positive for positive control) cells to wells in bottom of receiver plate (10,000 cells/well)
- Carefully place insert on top of the receiver plate
- Add 35 µl of BiTE® dilution to top chamber of insert
- Add 40 µl target-positive cells/T cell mixture to top chamber of insert
- Allow plates to equilibrate for 30 min. at room temperature
- Incubate at 37°C, 5% CO₂ for 72 hours

3. Perform cytotoxicity assay on the cells in the bottom chambers/receiver plates (imaging or CellTiter-Glo® - see TDCC protocol).