

# Propagating T5-phages for Fluorescent Staining

Dr. Steven Wilhelm

## Abstract

Please contact Dr. Steven Wilhelm (wilhelm@utk.edu) for additional information regarding this protocol.

**Citation:** Dr. Steven Wilhelm Propagating T5-phages for Fluorescent Staining. **protocols.io**

dx.doi.org/10.17504/protocols.io.igzcbx6

**Published:** 24 Jun 2017

## Protocol

### Propagating T5-phages

#### Step 1.

Grow a culture of E.coli ATCC11303 in LB media overnight

### Propagating T5-phages

#### Step 2.

The next morning, prewarm LB-agar plates to 37°C

### Propagating T5-phages

#### Step 3.

Pellet 1 mL of the cells in an Eppendorf tube at 10,000 xg, 2 min

### Propagating T5-phages

#### Step 4.

Discard the medium

### Propagating T5-phages


#### Step 5.

Resuspend the pellet into 1 mL 10 mM MgSO<sub>4</sub>

#### AMOUNT

1 ml Additional info:

#### REAGENTS

 Magnesium sulfate heptahydrate by Contributed by users

### Propagating T5-phages

## Step 6.

Melt 0.6% top agar in the microwave

Propagating T5-phages

## Step 7.

Prepare sterile 5 mL tubes and pipet 3 mL of melted top agar into them

📄 AMOUNT

3 ml Additional info:

Propagating T5-phages

## Step 8.

Keep tubes in 45°C heat block or water bath so that the agar stays in a liquid state, but not too hot to kill the cells.

Propagating T5-phages

## Step 9.

Add 10 µL T5 from the stock (except the control) to 100 µL of resuspended cells

Propagating T5-phages

## Step 10.

Incubate 5 min at room temperature

Propagating T5-phages

## Step 11.

Add the infected cells to the molton top agar

Propagating T5-phages

## Step 12.

Vortex briefly

Propagating T5-phages

## Step 13.

Pour onto prewarmed LB agar plate and tilt the plate to spread

Propagating T5-phages

## Step 14.

Let top agar solidify 30-60 min

Propagating T5-phages

## Step 15.

Evert plates and return to 37°C incubator

#### Propagating T5-phages

##### Step 16.

Let grow 6-8 hrs

#### Propagating T5-phages

##### Step 17.

When confluent plaques are seen on the plate, add 5 mL of sterile phage buffer (10 mM MgSO<sub>4</sub>; 1 mM CaCl<sub>2</sub>; 10 mM Tris-HCl, pH 7.5; 1% gelatin)

#### Propagating T5-phages

##### Step 18.

Let the plates sit overnight at 4°C with the buffer on top to gather the phages

#### Propagating T5-phages

##### Step 19.

The next day, transfer the liquid from the plates to a sterile tube and add 100 µL chloroform and 5 mL phage buffer

#### Propagating T5-phages

##### Step 20.

Shake for 1 min to let chloroform settle out

#### Propagating T5-phages

##### Step 21.

Transfer 4 mL of the supernatant to a fresh tube

#### AMOUNT

4 ml Additional info:

#### Propagating T5-phages

##### Step 22.

Filter the solution with 0.22 µm pore size syringe filter to get rid of all the cellular debris

#### Propagating T5-phages

##### Step 23.

Store phage at 4°C. They will stay infective for >1 yr.

#### PROTOCOL

## . [Fluorescent Staining of T5-phages](#)

CONTACT: [Steven Wilhelm](#)

Preparing T5 for Staining

### Step 23.1.

Pellet phage in an ultracentrifuge at 28,000 rpm, 90 min, 10°C in 35 mL polycarbonate tubes, filled to the rim, balanced with phage buffer so that the difference is <0.01 g.

#### NOTES

**Alyssa Alsante** 23 Jun 2017

Use the SW-28 rotor, or equivalent, cooled down in the fridge prior to using.

Preparing T5 for Staining

### Step 23.2.

Discard the supernatant

Preparing T5 for Staining

### Step 23.3.

Resuspend the pellet with phage buffer

#### AMOUNT

800 µl Additional info:

Preparing T5 for Staining

### Step 23.4.

If the pellet still seems to have some cellular debris in it, filter the resuspension using a 1 mL syringe and a small 0.2 µm pore size syringe filter.

Preparing T5 for Staining

### Step 23.5.

Re-pellet the viruses in a Beckmann TL-100 ultracentrifuge at 35,000 rpm, 2 hr, 6°C using the TL-55 rotor cooled down before use using 1.4 mL polycarbonate tubes, filled and balanced as the bigger tubes in step 1.

Preparing T5 for Staining

### Step 23.6.

Resuspend the pellet with phage buffer and transfer into 1.5 mL screw cap tubes

#### AMOUNT

100 µl Additional info:

Preparing T5 for Staining

### Step 23.7.

Store at 4°C

#### Staining T5 for Staining

##### Step 23.8.

Add YO-PRO stain to each tube containing 100 µL resuspension

##### AMOUNT

0.5 µl Additional info:

##### NOTES

**Alyssa Alsante** 23 Jun 2017

Make sure to work in the dark, because the stain is photosensitive.

#### Staining T5 for Staining

##### Step 23.9.

Cover the tubes with aluminum foil and let sit at 4°C for 48 hrs

#### Staining T5 for Staining

##### Step 23.10.

After incubation, increase the volume to 1.4 mL

#### Staining T5 for Staining

##### Step 23.11.

Centrifuge in the Beckmann TL-100 as in step 5

#### Staining T5 for Staining

##### Step 23.12.

Discard the supernatant

#### Staining T5 for Staining

##### Step 23.13.

Suspend the viruses in Milli-Q water

##### AMOUNT

100 µl Additional info: