



Feb 01, 2020

Lipid Annotation of MALDI IMS Datasets

Elizabeth Neumann¹, Jamie Allen¹, Jeff Spraggins¹, Danielle Gutierrez¹

¹Vanderbilt University

1 Works for me dx.doi.org/10.17504/protocols.io.864hzgw

VU Biomolecular Multimodal Imaging Center

Human BioMolecular Atlas Program (HuBMAP) Method Development Community



Jamie Allen ⚡

ABSTRACT

Scope:

Annotate lipid species detected by MALDI IMS analysis.

Expected Outcome:

A list of lipid identifications with low mass errors (<5 ppm) that correlate to ion images produced by MALDI IMS analysis.

- 1 Perform a quadratic recalibration of the IMS dataset using common well-characterized lipids:
 m/z 734.5694, 760.5851, 782.5670, 788.6164, 798.5410, 810.6007, 832.5827, and 848.5566
Software used are DataAnalysis or Mmass, and ppm error calculated in excel.
- 2 Create mass list from the averaged, re-calibrated IMS data set.
- 3 Calculate ppm error associated with common lipids, such as $[PC(32:0)+H]^+$ and $[PC(34:1)+H]^+$, for mass errors. If mass error is above 3 ppm, re-calibrate.
- 4 If LC-MS/MS has been performed, annotate each m/z value against the LC-MS/MS list, otherwise move to 4.1
- 4.1 Use LIPIDMAPS database to annotate each m/z value with the following search criteria:
https://www.lipidmaps.org/resources/tools/bulk_structure_searches.php?database=COMP_DB
 - 1a. Positive Mode Adducts: $[M+H]^+$, $[M+H-H_2O]^+$, $[M+Na]^+$, $[M+K]^+$
 - 1b. Negative Mode Adducts: $[M+H]^-$, $[M+Cl]^-$, $[M+HCOO]^-$, $[M+OAc]^-$
 2. Specify Mass Tolerance: 0.005 m/z
 3. Specify Chains: Even Chains Only
 4. Sort: Delta
- 5 Calculate ppm error associated with each assignment and remove assignments with errors larger than 5 ppm, although most have errors lower than 3 ppm.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited