

# Transformation of *Synechococcus* sp PCC7002

Anne Vogel, Rahmi Lale, Martin Hohmann-Marriott

## Abstract

This protocol describes the transformation of *Synechococcus* sp. PCC7002 with linear DNA fragments with homologous regions.

**Citation:** Anne Vogel, Rahmi Lale, Martin Hohmann-Marriott Transformation of *Synechococcus* sp PCC7002. **protocols.io** dx.doi.org/10.17504/protocols.io.kg2ctye

**Published:** 16 Jul 2018

## Guidelines

Work as sterile as possible.

## Protocol

### Culture Preparation

#### Step 1.

- Set up a starter culture of *Synechococcus* (e.g. from a -80°C stock or from a plate or a liquid culture) in AA+ medium (33°C incubation with bubbling, ca. 150uE)
- Incubate until OD<sub>730</sub> is above 1

#### Step 2.

- Use the starter culture to set up fresh culture in the bubbling flasks (cleaned and sterile) with an OD<sub>730</sub> of 0.07 in 150mL AA+ medium
- Incubate at 33°C with bubbling until OD<sub>730</sub> of 0.4 (ca. 18 hours in our lab)

#### Step 3.

- Transfer cells from the bubbling flasks into two sterile 50mL tubes
- Spin down the cells at 2500g (4700rpm) for 8 minutes at room temperature (22°C)
- Carefully discard the AA+ medium

### Transformation

#### Step 4.

- Resuspend cells in AA+ medium to a final OD<sub>730</sub> of 8
- Place 100 of the resuspended cells in sterile test tubes (15mL volume)
- Add 1µg of linear DNA with 50-500nt homologous regions (e.g. your purified PCR product)

#### Step 5.

- Place tubes in a rack in a 30°C cabinet with low light (ca. 30uE) for six hours
- In the meantime, place sterile filters (Nucleopore SN 145318) on **non-selective** AA+ plates with sterile tweezers.

**Step 6.**

- Plate the whole volume of cells onto the sterile filters on the AA+ plates
  - Spread the cells carefully
  - Leave the plates to dry in the flow cabinet (30 minutes)
- 

**Step 7.**

- Grow the cells overnight in the 30°C cabinet with low light (ca. 30 uE)
- 

**Step 8.**

- Transfer the filters with the cells onto AA+ plates with the appropriate antibiotics
  - Grow cells in the 30°C cabinet with low light (ca. 30uE) for 2 days.
- 

**Step 9.**

- Transfer cells to 33°C incubator with higher light (ca. 150uE)
  - Place the plates near the door so that they do not get too much light.
  - Be careful to seal the plates well with parafilm so that they do not dry out.
- 

**Step 10.**

- Wait until colonies appear on the filter. A lot of cells will die, so there is most likely a yellowish background on which the green colonies will appear.
- 

## Segregation

**Step 11.**

- Pick single colonies with a sterile pipettip and streak them onto selective AA+ plates.
  - Incubate plates in the 33C cabinet the same way as done before with the selective plates.
- 

**Step 12.**

- Wait until colonies appear and repeat the procedure
  - To be sure that complete segregation occurred, perform colony PCR on your colonies.
-