

# Increased sensitivity of Euplotes crassus to selective agents using 0.3 M glucose-based culture conditions.

#### Lawrence A. Klobutcher

#### **Abstract**

Numerous unpublished studies indicate that *Euplotes crassus* is resistant to extremely high concentrations of commonly used selective agents when grow in artificial seawater. We have found that when *E. crassus* is grown in 1 part artificial seawater + 9 parts 0.3 M glucose, they are much more sensitive to a number of selective agents. For example, 1080 ug/ml paromomycin in artificial seawater only inhibited growth, as little as 120 ug/ml was effective in killing cells when grown with 0.3 M glucose. Similarly, G418 at 400 ug/ml had little effect in seawater, but as little as 100 ug/ml blocked growth in 0.3 M glucose. This may not prove true for all selective agents, as we saw no effect with paclitaxel up to 100 uM under either growth condition.

**Citation:** Lawrence A. Klobutcher Increased sensitivity of Euplotes crassus to selective agents using 0.3 M glucose-based culture conditions.. **protocols.io** 

dx.doi.org/10.17504/protocols.io.gwnbxde

Published: 11 Jan 2017

#### **Protocol**

#### Step 1.

Distribute a dense culture of *Dunaliella salina* to two 50 ml screw-cap centrifuge tubes, and pellet the algae by centrifugation at 1,000 RPM (200 x g) for 2 min. in a clinical centrifuge.

#### Step 2.

Remove all but the last 2.5 ml of supernatant from each tube, resuspend the algae pellets, and combine.

#### Step 3.

To the above tube, add 2.5 ml of a log phase culture of *Euplotes crassus* (grown with *Dunaliella* in artificial seawater as the food source), and 45 ml of 0.3 M glucose (freshly prepared using deionized water).

#### Step 4.

Add selective agent to desired final concentration. 200 ug/ul G418, 200 ug/ml paromomycin, or 40 ug/ml puromycin.

#### Step 5.

Cells can be distributed microtiter plates. For example, 2 ml per well for a 24-well microtiter plate.

## Step 6.

Incubate at room temperature in a humidified chamber (e.g., plastic box with wet paper towels or beaker of water). Cell death occurs within 2-3 days.

### Step 7.