

# Euplotes crassus transformation using Lipofectamine 2000 as vehicle Version 2

RACHELE CESARONI

## Abstract

**Citation:** RACHELE CESARONI Euplotes crassus transformation using Lipofectamine 2000 as vehicle. **protocols.io**  
dx.doi.org/10.17504/protocols.io.g3xbypn

**Published:** 30 Jan 2017

## Protocol

### Step 1.

$2 \times 10^4$  well-fed cells (we used E.coli as the only food source) were collected by centrifugation (400 rcf for 3 minutes).

### Step 2.

Cells were washed twice with artificial sea water (see attachment for the recipe) and once with 500 mM sorbitol, 0.5 mM Tris-HCl pH 7.0 (400 rcf for 3 minutes each time), in which they were also resuspended in 50  $\mu$ l (to get this small volume an additional minute of centrifugation may be required).

## ⊕ NOTES

**Estienne Swart** 30 Jan 2017

Recipe for complete sea water (1 L):

36 g Reef Crystals

1 ml Walne's solution

1 ml of 10  $\mu$ g/ml  $\text{FeSO}_4$

0.2 ml of 2 mg/ml thiamine (light sensitive; store at 4°C)

Add distilled water to 1 L

### Step 3.

Cells were transferred into two wells within 96-well plates for transfection.

## ■ ANNOTATIONS

**Rachele Cesaroni** 30 Jan 2017

One of the two wells was used for the negative control.

### **Step 4.**

2.5 µl of Lipofectamine 2000 Reagent were diluted in 25 µl of the same medium of the cells (500 mM sorbitol, 0.5 mM Tris-HCl pH 7.0).

### **Step 5.**

5 µg of DNA (0.5-5 µg/µl) dissolved in MilliQ H<sub>2</sub>O were diluted in 125 µl of the same medium of the cells (500 mM sorbitol, 0.5 mM Tris-HCl pH 7.0).

### **Step 6.**

25 µl of the diluted DNA were added to 25 µl of the diluted Lipofectamine 2000 Reagent (1:1 ratio) and incubated for 10 minutes at room temperature (before incubation mix by pipetting up and down 5 times or vortex for 10 seconds).

### **Step 7.**

10 µl of the transfection complexes were added drop-wise to the cells in the well plates (gently swirl the dish to ensure uniform distribution of the transfection complexes).

### **Step 8.**

One hour after addition of Lipofectamine 2000 complexes, 50 µl of artificial sea water were added to the cells.

### **Step 9.**

After another hour other 50 µl of artificial sea water were added to the cells.

### **Step 10.**

An hour later cells were harvested (400 rcf for 3 minutes) and washed twice with artificial sea water (400 rcf for 3 minutes each time), in which they were then also resuspended in 400 µl.

### **Step 11.**

Cells were transferred into depression wells.

### **Step 12.**

Cells were incubated at least for three days at 24°C, then analyzed by fluorescence microscopy to determine gene expression.