

Protein extraction, alkylation, and digestion for LC/MS of HEK-293 Version 2

Shun Adachi

Abstract

Protocol for protein extraction, alkylation, and digestion to obtain LC/MS (liquid chromatography-mass spectrometry) data of calculated protein signal intensities in HEK-293 cells to monitor them.

Citation: Shun Adachi Protein extraction, alkylation, and digestion for LC/MS of HEK-293. **protocols.io**
dx.doi.org/10.17504/protocols.io.h59b896

Published: 25 May 2017

Before start

Prepare a Krebs-Ringer-Buffer:

- 154 mM NaCl
- 5.6 mM KCl
- 5.5 mM glucose
- 20.1 mM HEPES pH 7.4
- 25 mM NaHCO₃

Materials

- RIPA buffer [08714](#) by [Nacalai Tesque](#)
- ✓ NaCl by Contributed by users
- ✓ KCl by Contributed by users
- ✓ Glucose by Contributed by users
- ✓ HEPES pH 7.4 by Contributed by users
- ✓ NaHCO₃ by Contributed by users
- Micro BCA Protein Assay Kit [23235](#) by [Thermo Fisher Scientific](#)
- XL-Tryp Kit Direct Digestion [YC - 3120](#) by [APRO SCIENCE](#)
- In-Gel R-CAM Kit [YC-5004](#) by [APRO SCIENCE](#)
- ZipTipC18 by [Merck Millipore](#)

Protocol

Protein extraction

Step 1.

Extract the HEK-293 proteins using the standard protocol for the [RIPA buffer \(NACALAI, INC., Kyoto, Japan\)](#).

Alkylation and Digestion

Step 2.

Wash approximately 10^6 harvested cells once in Krebs-Ringer-Buffer (KRB; 154 mM NaCl, 5.6 mM KCl, 5.5 mM glucose, 20.1 mM HEPES pH 7.4, 25 mM NaHCO_3).

Alkylation and Digestion

Step 3.

Resuspend the cells in 30 μl of RIPA buffer, passed in and out through 21G needles for destruction, and incubate on ice for 1 h.

 DURATION

01:00:00

Alkylation and Digestion

Step 4.

Centrifuge the cells at 10,000 g for 10 min at 4°C.

 DURATION

00:10:00

Alkylation and Digestion

Step 5.

Collect the supernatants.

Alkylation and Digestion

Step 6.

Quantify the proteins by using a Micro BCA Protein Assay Kit (Thermo Fisher Scientific, Waltham, U.S.A.).

Alkylation and Digestion

Step 7.

Solidify the samples in acrylamide gel. For further steps, use XL-Tryp Kit Direct Digestion (APRO SCIENCE, Naruto, Japan).

Alkylation and Digestion

Step 8.

Wash the samples in ultrapure water. (1/2)

Alkylation and Digestion

Step 9.

Wash the samples in ultrapure water. (2/2)

Alkylation and Digestion

Step 10.

Wash the samples in dehydration solution. (1/3)

Alkylation and Digestion

Step 11.

Wash the samples in dehydration solution. (2/3)

Alkylation and Digestion

Step 12.

Wash the samples in dehydration solution. (3/3)

Alkylation and Digestion

Step 13.

Dry the samples.

Alkylation and Digestion

Step 14.

Process the samples using an In-Gel R-CAM Kit (APRO SCIENCE, Naruto, Japan).

Alkylation and Digestion

Step 15.

Reduce the samples for 2 h at 37°C.

 DURATION

02:00:00

Alkylation and Digestion

Step 16.

Alkylate the samples for 30 min at room temperature.

 DURATION

00:30:00

Alkylation and Digestion

Step 17.

Wash samples with ultrapure water. (1/5)

Alkylation and Digestion**Step 18.**

Wash samples with ultrapure water. (2/5)

Alkylation and Digestion**Step 19.**

Wash samples with ultrapure water. (3/5)

Alkylation and Digestion**Step 20.**

Wash samples with ultrapure water. (4/5)

Alkylation and Digestion**Step 21.**

Wash samples with ultrapure water. (5/5)

Alkylation and Digestion**Step 22.**

Wash samples with destaining solution. (1/2)

Alkylation and Digestion**Step 23.**

Wash samples with destaining solution. (2/2)

Alkylation and Digestion**Step 24.**

Dry the samples.

Alkylation and Digestion**Step 25.**

Trypsinize the resultant samples overnight at 35°C.

Alkylation and Digestion-Day2**Step 26.**

Collect the dissolved digested peptides by ZipTipC18 (Merck Millipore Corp., Billerica, U.S.A.).

Alkylation and Digestion-Day2

Step 27.

Dampen the tips with acetonitrile x three tips. (1/2)

Alkylation and Digestion-Day2

Step 28.

Dampen the tips with acetonitrile x three tips. (2/2)

Alkylation and Digestion-Day2

Step 29.

Equilibrate tips with 0.1% trifluoroacetic acid x three tips. (1/2)

Alkylation and Digestion-Day2

Step 30.

Equilibrate tips with 0.1% trifluoroacetic acid x three tips. (2/2)

Alkylation and Digestion-Day2

Step 31.

Collect the peptides by 20 cycles of aspiration and dispensing x three tips.

Alkylation and Digestion-Day2

Step 32.

Wash the peptides with 0.1% trifluoroacetic acid x three tips. (1/2)

Alkylation and Digestion-Day2

Step 33.

Wash the peptides with 0.1% trifluoroacetic acid x three tips. (2/2)

Alkylation and Digestion-Day2

Step 34.

Elute by 0.1% trifluoroacetic acid /50% acetonitrile with aspiration and dispense x three tips. (1/5)

Alkylation and Digestion-Day2

Step 35.

Elute by 0.1% trifluoroacetic acid /50% acetonitrile with aspiration and dispense x three tips. (2/5)

Alkylation and Digestion-Day2

Step 36.

Elute by 0.1% trifluoroacetic acid /50% acetonitrile with aspiration and dispense x three tips. (3/5)

Alkylation and Digestion-Day2

Step 37.

Elute by 0.1% trifluoroacetic acid /50% acetonitrile with aspiration and dispense x three tips. (4/5)

Alkylation and Digestion-Day2

Step 38.

Elute by 0.1% trifluoroacetic acid /50% acetonitrile with aspiration and dispense x three tips. (5/5)

Alkylation and Digestion-Day2

Step 39.

Vacuum dry the peptides.

Alkylation and Digestion-Day2

Step 40.

Store the finalized samples at -20°C .

LC/MS

Step 41.

Before performing LC/MS, resuspend samples in 0.1% formic acid.

LC/MS

Step 42.

Quantify the amounts with the Pierce Quantitative Colorimetric Peptide Assay (Thermo Fisher Scientific, Waltham, U.S.A.).