



Dec 06, 2019

## Dye-terminator DNA sequencing

[Diep R. Ganguly<sup>1</sup>](#)<sup>1</sup>The Australian National University

1

Works for me

[dx.doi.org/10.17504/protocols.io.97jh9kn](https://dx.doi.org/10.17504/protocols.io.97jh9kn)[Diep R. Ganguly](#)

The Australian National University



### ABSTRACT

This protocol (based on the BigDye® Terminator v3.1 Cycle Sequencing Kit) is for performing terminator cycling sequencing reactions for Sanger sequencing of amplified PCR products or plasmid DNA on the 3130X genetic analyser (Applied Biosystems).

□  
BigDye Terminator  
v3.1.pdf

### MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Antarctic Phosphatase - 1,000 units	M0289S	New England Biolabs
96 well PCR Plate Non-skirted	MPS-499	Phenix Research
Nuclease-free water (e.g. MilliQ or HPLC grade water)		
primers		
EDTA		
10 mM dNTPs	10297-018	Life Technologies
Ethanol	100983	Merck Millipore
BigDye™ Terminator v3.1 Cycle Sequencing Kit	4337454	Thermo Fisher
Exonuclease I (E. coli)	M0293S	NEB
Hi-Di™ Formamide	4311320	Thermo Fisher Scientific

### BEFORE STARTING

Optimize PCR cycling (if sequencing amplified PCR products) to ensure your reaction produces a single product. Perform gel excision or PCR clean-up with the potential inclusion of incubating with Antarctic phosphatase and Exonuclease 1 to dephosphorylate and degrade unincorporated dNTPs in PCR reactions to prepare templates for DNA sequencing.

### Terminator cycling reaction

- 1 Perform sequencing reaction with BigDye Terminator cycling kit and either forward or reverse primers, or both.

Component	Volume (µl)
2.5X Reaction Ready Mix	4
5X BigDye Sequencing buffer	2
20 µM F/R Primer	1
Template (plasmid or cleaned PCR product)	150-300 ng dsDNA or approx 10ng PCR product (see BigDye manual)
Nuclease-free water	to 20 µl

### BigDye Terminator Cycling reaction

- 2 Set the following protocol and allow thermal cycler to reach 96 °C:
  1. 1 min at 96 °C
  2. 30 cycles consisting of: 96 °C for 10 seconds, 50 °C for 5 seconds, and 60 °C for 4 min.
  3. Hold at 10 °C.

### Purify products

- 3 To the PCR reaction, add 60 µL 100% ethanol and 5 µl 125 mM EDTA.
- 4 Incubate at room temperature for 15 minutes.
- 5 Centrifuge at 4 °C at max speed for 10 minutes.
- 6 Remove as much supernatant as possible, then allow to air-dry in the dark for 10-15 minutes.

### Resuspend products and submit for sequencing

- 7 Resuspend the pellet (likely transparent) in 20 µL HiDi Formamide (add 20 µL to any empty wells). Spin down briefly.
- 8 Submit for sequencing on 3130X genetic analyser (Applied Biosystems).



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited