



Jul 11, 2019

## Heterologous expression and affinity purification of Strep-tagged (KaiC) proteins

Anika Wiegard<sup>1</sup>, Christin Köbler<sup>2</sup>, Katsuaki Oyama<sup>3</sup>, Anja K. Dörrich<sup>4</sup>, Chihiro Azai<sup>5,3</sup>, Kazuki Terauchi<sup>5,3</sup>, Annegret Wilde<sup>2</sup>, Ilka Maria Axmann<sup>6</sup>

<sup>1</sup>Karolinska Institutet, Department of Cell and Molecular Biology, Stockholm, Sweden, <sup>2</sup>Institute of Biology III, Faculty of Biology, University of Freiburg, 79104 Freiburg, Germany, <sup>3</sup>Graduate School of Life Sciences, Ritsumeikan University, Kusatsu, Shiga 525-8577, Japan, <sup>4</sup>Institute for Microbiology and Molecular Biology, Justus-Liebig University, 35392 Giessen, Germany, <sup>5</sup>College of Life Sciences, Ritsumeikan University, Kusatsu, Shiga 525-8577, Japan, <sup>6</sup>Institute for Synthetic Microbiology, Cluster of Excellence on Plant Sciences (CEPLAS), Heinrich Heine University Duesseldorf, 40225 Duesseldorf, Germany

1

Works for me

dx.doi.org/10.17504/protocols.io.meac3ae

Axmann Lab

CyanoWorld



Anika Wiegard

Karolinska Institutet, Department of Cell and Molecular Biol...



### ABSTRACT

This protocol can be used to:

- (i) express Strep-tagged proteins in *E. coli*
- (ii) lyse cells
- (iii) purify Strep-tagged proteins via gravity flow affinity chromatography using either Strep-tactin XT superflow or Strep-Tactin resin

### MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
NaOH		
Magnesium chloride hexahydrate	<a href="#">View</a>	<a href="#">Sigma Aldrich</a>
Tris(hydroxymethyl)aminomethane	252859-500G	<a href="#">Sigma Aldrich</a>
NaCl		
Lysozyme	12671-19-1	<a href="#">Sigma Aldrich</a>
Benzonase	101654	<a href="#">Merck Millipore</a>
Roche Complete Protease Inhibitor EDTA-Free tablets	5056489001	<a href="#">Sigma Aldrich</a>
Ampicillin sodium salt	A0166	<a href="#">Millipore Sigma</a>
HCl	<a href="#">View</a>	
14 Dithiotreitol (DTT)	6908.1	<a href="#">Carl Roth</a>
Adenosin-5-triphosphate disodium salt (ATP)	HN35.1	<a href="#">Carl Roth</a>
Quick Start™ Bradford 1x Dye Reagent	5000205	<a href="#">Bio-rad Laboratories</a>
LB broth (Lennox)	X964.1	<a href="#">Carl Roth</a>
pASK-IBA5plus	2-1404-000	
Rosetta-gami™ B (DE3) Competent Cells - Novagen	71136	
Chloramphenicol	3886.2	<a href="#">Carl Roth</a>
D( )Biotin	2-1016-002	

NAME ▾	CATALOG # ▾	VENDOR ▾
Anhydrotetracycline	2-0401-001	

  

NAME ▾	CATALOG # ▾	VENDOR ▾
Strep-Tactin®XT Superflow® 50% suspension	2-4010-010	
Strep-Tactin® Sepharose® 50% suspension	2-1201-010	
D( )Biotin	2-1016-002	
Desthiobiotin	2-1000-001	
10x Buffer R Strep-Tactin® Regeneration Buffer with HABA	2-1002-100	

## heterologous protein expression in E.coli

- 1 transformation
  - transform *E.coli* expression cells (e.g. Rosetta gamiB (DE3), or Rosetta gami2 (DE3)) with your pASK based expression plasmid
- 2 pre-culture:
  - inoculate LB medium containing 75-100 µg ampicillin ml<sup>-1</sup> (optional: plus 30 µg chloramphenicol ml<sup>-1</sup>) with resulting transformants (use e.g. 200 ml LB)
  - incubate over night at 37 °C and 200-250 r.p.m.
- 3 expression culture:
  - inoculate LB medium containing 75-100 µg ampicillin ml<sup>-1</sup> (optional: plus 30 µg chloramphenicol ml<sup>-1</sup>) with 4-10 % volume of the pre-culture (add e.g. 200 ml pre-culture to 1.8 l LB)
  - incubate at 37 °C and 200-250 r.p.m. until OD<sub>600nm</sub> = ~0.5

*Note: use erlenmeyer flasks with a volume of at least 4l to ensure sufficient aeration.*
- 4 induction of protein expression:
  - when cell density reaches OD<sub>600nm</sub> = ~0.3-0.7 (optimal 0.5) add 200 µl of 2 mg anhydrotetracycline ml<sup>-1</sup> (final concentration 200 ng anhydrotetracycline ml<sup>-1</sup>)
  - choose optimal expression condition (has to be tested for each protein of interest), e.g.: over night at 18-25 °C and 200-250 r.p.m. or at 35-37 °C and 200-250 r.p.m. for 3.5-5.5 hours

## cell harvest

- 5
  - spin cells for 10 min at 4°C and 4000g
  - discard supernatant
  - keep cells on ice

🕒 00:10:00 centrifugation

## cell disruption

- 6 enzymatic lysis by lysozyme:
  - weight cell pellets
  - resuspend cells in ice-cold buffer W [20mM Tris/HCl (pH8), 150 mM NaCl, 2 mM DTT (*only for KaiC proteins: 5 mM MgCl<sub>2</sub>, 1 mM ATP*)] including protease inhibitors (e.g. protease inhibitor cocktail, Roche) using a paint brush. (use 3 ml buffer per g cells)
  - add a spatula tip's worth of lysozyme (or add lysozyme stock solution to a final concentration of 1 mg lysozyme ml<sup>-1</sup>)
  - add 125 U benzonase
  - incubate on ice for 30 min

*Note: you can also use 50 mM Tris in buffer*

Note: you can skip enzymatic lysis step and perform longer sonication instead

🕒 00:30:00 incubation on ice

## 7 sonication:

- sonicate the cell suspension for 6 min on ice using e.g. a Bandelin sonopuls homogenizer and the following parameters
- Note: During sonication, the temperature of the cell suspension should be kept below 15 °C

tip	KE76
cycle	1 (0.1 sec active cycle, 0.9 sec passive cycle)
output	60 %

Note: if you did not perform enzymatic lysis as described before, use alternative prolonged sonication conditions (e.g. alternation of 10 sec pulse and 10 sec pause for 25 minutes with 30 % output)

🕒 00:06:00 sonication

## 8 clarification of the lysate:

- centrifuge the resulting lysate for 20 min at 4 °C and 23000 g to remove insolubles
- keep the resulting supernatant (= soluble proteins)

Note: Many conical centrifugation tubes cannot withstand centrifugation of 23000g. If you want to use them, you can reduce centrifugal force, while increasing centrifugation time.

🕒 00:20:00 centrifugation

## affinity purification

## 9 preparation of affinity column (at 4°C or room temperature)

- pour 3 ml Strep-tactin XT superflow resin (or Strep-Tactin Sepharose) in an appropriate glass column
- equilibrate with 30 ml ice-cold buffer W
- alternatively: equilibrate with 25 ml ice-cold buffer W and subsequently with 5 ml ice-cold buffer W including protease inhibitor

Note: if resin was stored in buffer R before, the colour will change from orange to white



Strep-Tactin®XT Superflow® 50%  
suspension  
Catalog #: 2-4010-010



Strep-Tactin® Sepharose® 50% suspension  
Catalog #: 2-1201-010

## 10 protein binding (at 4°C or room temperature)

- apply soluble proteins to your column
- collect the flow through

## 11 washing (at 4°C or room temperature)

- wash column with 15-50 ml ice-cold buffer W including protease inhibitors (at 4°C or room temperature)

## 12 elution from Strep-Tactin XT superflow resin (if you use Strep-Tactin Sepharose proceed with step 13 instead) (at 4°C or room temperature)

- elute proteins with 9 ml ice cold buffer BXT (buffer W + 50 mM D(+)biotin)
- collect 6 fractions of 1.5 to 2ml



D(+)Biotin  
Catalog #: 2-1016-002

13 elution from Strep-Tactin Sepharose (if you use Strep-tactin XT superflow resin skip this step and move on to step 14) (at 4°C or room temperature)

- elute proteins with 15-30 ml ice cold buffer W + 2.5 mM desthiobiotin
- collect fractions of 1 to 2ml



Desthiobiotin

Catalog #: 2-1000-001

14 regeneration of Strep-Tactin XT superflow resin (if you use Strep-Tactin Sepharose proceed with step 15 instead) (at room temperature)

- wash with 15 ml 10 mM freshly prepared NaOH
- optional: regeneration can be tested by adding buffer R (100 mM Tris-Cl, 150 mM NaCl, 1 mM EDTA, 1 mM HABA (hydroxy-azophenyl-benzoic acid)). An orange-shift indicates successful regeneration

*Note: you cannot use buffer R for regeneration of Strep-tactin XT superflow resin*



10x Buffer R Strep-Tactin® Regeneration

Buffer with HABA

Catalog #: 2-1002-100

15 regeneration of Strep-Tactin Sepharose (if you use Strep-tactin XT superflow resin, you already finished regeneration with step 14) (at room temperature)

- wash with buffer R (0.1 M Tris-Cl, 150 mM NaCl, 1 mM EDTA, 1 mM HABA (hydroxy-azophenyl-benzoic acid) until agarose turns orange

16 qualitative analysis of eluate fractions:

- for each fraction, pipette 80 µl of Bradford solution in a well of a 96 well plate
- add 5-20 µl of your fraction
- colour change to blue indicates successful elution of proteins --> keep those fractions
- control quality of eluted protein by separation via SDS-PAGE (optional: analyse aliquots of lysate, flow through and washing steps in parallel)

17 buffer exchange:

- a) perform size exclusion chromatography or
- b) exchange buffer using centrifugal concentrators
  - mix all eluate fractions of sufficient protein quality and transfer them to a disposable centrifugal concentrator
  - concentrate protein by centrifugation
  - add your desired buffer
  - concentrate again
  - repeat this step until the buffer is completely exchanged

*Note: choose the molecular cut-off of the concentrator and centrifugal force according to the manufacturer's instructions.*



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited