



May 10,
2019

Working

U Mass - Cholesterol (HDL) [↗](#)

Jason Kim¹

¹University of Massachusetts

dx.doi.org/10.17504/protocols.io.xt6fnre

Mouse Metabolic Phenotyping Centers

Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

Summary:

This experiment involves a spectrophotometric measurement using Roche Cobas Clinical Chemistry Analyzer. Serum levels of HDL-cholesterol reflect cholesterol metabolism and are associated with cardiovascular disease.

EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=170&docType=Protocol>

MATERIALS

NAME	CATALOG #	VENDOR	CAS NUMBER	RRID
HDL Cholesterol Plus 3rd gen	05401488 190	Roche		
Calibrator f.a.s. Lipids	12172623 160	Roche		
Precinorm L	10781827 122	Roche		
Precipath HDL/LDL-C	11778552 122	Roche		
NaCl Diluent 9%	04774230 190	Roche		
Cleaner	04774248 190	Roche		
Micro Sample cups	11406680 001	Roche		
NERL High Quality Water	9805	Fisher Scientific		

MATERIALS TEXT

Note:

Roche, [RRID:SCR_001326](#)

Fisher Scientific, [RRID:SCR_008452](#)

BEFORE STARTING

Notes:

✓ Try to use freshly prepared serum and plasma samples for this assay.

✓ No dilution or treatment of the sample is required, but plasma samples should be centrifuged to remove any fibrin/fibrinogen clumps.

✓ Samples should be stored at 2-8°C for 24 hours prior to analysis. For longer periods, store samples at -70°C, and avoid repeated

freeze/thaw cycles.

✓ A 50 µl dead volume is required in addition to sample volume for multi-protein analysis (typically 1-5 µl).

- 1 Perform daily quality control assessment of instrumentation before analysis.
- 2 Load each sample into a specialized micro-sample cup for the clinical chemistry analyzer.
- 3 Select Cholesterol (HDL) test on display and run the analysis.
- 4 Collect and analyze the data.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited