



Sep 13, 2019

## K/2 Ian / K-ET V.2

Ian Probert<sup>1</sup><sup>1</sup>CNRS-Sorbonne Université, Station Biologique, Place G. Tessier 29680 Roscoff FRANCE**1** Works for me dx.doi.org/10.17504/protocols.io.7bxhipn

Roscoff Culture Collection

Daniel Vaultot  
Station Biologique, Roscoff, France

## ABSTRACT

Modified from K medium by Ian Probert

## Reference

Keller, M.D., Selvin, R.C., Claus, W. &amp; Guillard, R.R.L. 1987. Media for the culture of oceanic ultraphytoplankton. J. Phycol. 23:633–8.

## BEFORE STARTING

Please refer to our general recommendations to grow cultures :

<https://www.protocols.io/private/A48906DC1374AD6281495CB86A8F092F>

- 1 ■ To 994 ml of seawater (pH 8.2, adjusted with NaOH) add:

Quantity	Compound	Stock solution (sterile)	Final conc. in K medium
0.5ml	NaNO <sub>3</sub>	48.9542g/litre H <sub>2</sub> O	288µM
0.5ml	NH <sub>4</sub> Cl *	0.535g/litre H <sub>2</sub> O	5µM
0.5ml	KH <sub>2</sub> PO <sub>4</sub>	4.8992g/litre H <sub>2</sub> O	18µM
0.5ml	FeEDTA solution	(see recipe below)	(see below)
0.5ml	Trace metal solution	(see recipe below)	(see below)
1.0ml	f/2 vitamin solution	(see recipe below)	(see below)

\* optional

## FeEDTA solution

- 2

- To 950ml distilled H<sub>2</sub>O add:

Quantity	Compound	Stock solution	Final conc. in K medium
4.3g	(Na)FeEDTA	-	5.85µM

- Make up to 1 litre with milliQ H<sub>2</sub>O, sterilize (filter 0.22µm) and store in fridge.

### Trace metal solution

- 3 ■ To 950ml distilled H<sub>2</sub>O add:

Quantity	Compound	Stock solution	Final conc. in K medium
37.22g	Na <sub>2</sub> EDTA.2H <sub>2</sub> O	-	50μM
1.0ml	CuSO <sub>4</sub> .5H <sub>2</sub> O	2.497g/litre H <sub>2</sub> O	0.005μM
1.0ml	Na <sub>2</sub> MoO <sub>4</sub> .2H <sub>2</sub> O	7.2585g/litre H <sub>2</sub> O	0.015μM
1.0ml	ZnSO <sub>4</sub> .7H <sub>2</sub> O	23.0g/litre H <sub>2</sub> O	0.004μM
1.0ml	CoSO <sub>4</sub> .7H <sub>2</sub> O	14.055g/litre H <sub>2</sub> O	0.025μM
1.0ml	MnCl <sub>2</sub> .4H <sub>2</sub> O	178.11g/litre H <sub>2</sub> O	0.45μM
1.0ml	H <sub>2</sub> SeO <sub>3</sub>	1.29g/litre H <sub>2</sub> O	0.005μM
1.0ml	NiCl <sub>2</sub> .6H <sub>2</sub> O	1.49 g/litre H <sub>2</sub> O	0.00314μM

- Make up to 1 litre with milliQ H<sub>2</sub>O, sterilize (filter 0.22μm) and store in fridge.

### f/2 Vitamin solution

- 4 ■ To 950ml distilled H<sub>2</sub>O add

Quantity	Compound	Stock solution	Final conc. in K medium
1.0ml	Vit. B <sub>12</sub> (cyanocobalamin)	0.5g/litre H <sub>2</sub> O	0.37nM
1.0ml	Biotin	5.0mg/litre H <sub>2</sub> O	2.0nM
100.0mg	Thiamine HCl	-	0.3μM



- Make up to 1 litre with milliQ H<sub>2</sub>O, filter sterilize into plastic vials and store in freezer.

### Sterilization of medium

- 5 ■ Optional: *Heat to 80°C for 2 hours and leave to cool – this should kill most organisms but should not chemically modify the medium too much*  
 ■ Filter sterilize through 0.22μm filters (e.g. Millipore Steritop units) into sterile (autoclaved) polycarbonate bottles.

### For K-ET

- 6 ■ Add 10-30 ml marine soil extract (ET) as detailed in the protocol below

Soil extract for algal media  
by Daniel Vaulot,  
Station Biologique, Roscoff, France

PREVIEW

RUN

- 6.1 Sample 10 g of dry soil from a place that does not contain any pesticide nor pollutant and where you can be sure to be able to come back in order to always use the same soil. The soil must be dry because the boiling step will be longer if the soil is not dry.
- 6.2 Add the 10 g soil to 400ml of milliQ water.
- 6.3 Boil during 1 hour.

6.4 Filter on 0.2 µm or first through 0.8/0.4 µm and then through 0.2µm in order not to clog the filter.

6.5 Aliquot and freeze.

6.6 Add 30 à 50ml of soil extract per 1 L of medium.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited