



Cantilever functionalisation for AFM single molecule force spectroscopy [↗](#)

PLOS One

Stefania Marcotti¹, Koichiro Maki², Gwendolen C Reilly¹, Damien Lacroix¹, Taiji Adachi³

¹University of Sheffield, ²University of Tokyo, ³Kyoto University

[dx.doi.org/10.17504/protocols.io.ra4d2gw](https://doi.org/10.17504/protocols.io.ra4d2gw)

Stefania Marcotti
University of Sheffield

ABSTRACT

Four different samples were employed for experiments, designed as follows:

1. HABP/HA: cantilever functionalised with hyaluronic acid binding protein, untreated cell sample;
2. BSA/HA: cantilever functionalised with bovine serum albumin, untreated cell sample;
3. untreated/HA: non-functionalised cantilever, untreated cell sample;
4. HABP/HAase: cantilever functionalised with hyaluronic acid binding protein, cell sample treated with HAase

The steps of activation are listed in this protocol and were the same for Sample 1 (HABP/HA, functionalisation molecule: hyaluronic acid binding protein), Sample 2 (BSA/HA, functionalisation molecule: bovine serum albumin) and Sample 4 (HABP/HAase, functionalisation molecule: hyaluronic acid binding protein). The cantilevers used to test cells in Sample 3 (untreated/HA) were not treated.

Similar protocols were described in:

* Maki K et al. Mechano-adaptive sensory mechanism of α -catenin under tension. Sci Rep 2016;6: 24878. doi:10.1038/srep24878

* Maki K et al. Real-time TIRF observation of vinculin recruitment to stretched α -catenin by AFM. Sci Rep 2018;8: 1–8.

doi:10.1038/s41598-018-20115-8

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0206056>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Marcotti S, Maki K, Reilly GC, Lacroix D, Adachi T (2018) Hyaluronic acid selective anchoring to the cytoskeleton: An atomic force microscopy study. PLoS ONE 13(10): e0206056. doi: [10.1371/journal.pone.0206056](https://doi.org/10.1371/journal.pone.0206056)

PROTOCOL STATUS







Working

MATERIALS

NAME	CATALOG #	VENDOR
MilliQ Water		
(3-Aminopropyl)triethoxysilane	440140	Sigma Aldrich
Bovine Serum Albumin	A9647	Sigma Aldrich
Maleimide-C3-NTA	M035	
2-Mercaptoethanol	M3148	Sigma Aldrich
Hyaluronic Acid Binding Protein, Bovine Nasal Cartilage, Biotinylated	385911	Merck Millipore
Trizma® hydrochloride solution	T2319	Sigma Aldrich


STEPS MATERIALS

NAME	CATALOG #	VENDOR
------	-----------	--------

NAME	CATALOG #	VENDOR
 (3-Aminopropyl)triethoxysilane	440140	Sigma Aldrich
MilliQ Water		
Maleimide-C3-NTA	M035	
 Trizma® hydrochloride solution	T2319	Sigma Aldrich
 Hyaluronic Acid Binding Protein, Bovine Nasal Cartilage, Biotinylated	385911	Merck Millipore
 Trizma® hydrochloride solution	T2319	Sigma Aldrich
MilliQ Water		
 Bovine Serum Albumin	A9647	Sigma Aldrich
 2-Mercaptoethanol	M3148	Sigma Aldrich
MilliQ Water		
MilliQ Water		

- 1 AFM cantilevers (OMCL-TR400PSA; spring constant, 0.02 N/m; curvature radius of tip, 15 nm; Olympus Co.) are oxidised using an ozone cleaner and submerged in 2% w/w (3-Aminopropyl)triethoxysilane / ultra-pure water for 15 minutes to depose (-SH) groups on the probe surface

• 2 Mass Percent (3-Aminopropyl)triethoxysilane / ultra-pure water

 (3-Aminopropyl)triethoxysilane
by Sigma Aldrich
Catalog #: 440140

 MilliQ Water


⌚ 00:15:00 (3-Aminopropyl)triethoxysilane

⚠ SAFETY INFORMATION


Please refer to SDS for (3-Aminopropyl)triethoxysilane

https://www.sigmaaldrich.com/catalog/product/aldrich/440140?lang=en&ion=GB&utm_medium=referral&utm_source=pubchem&utm_campaign=pubchem_2017

- 2 After washing, the cantilevers are submerged in 6 mM Maleimide-C3-NTA / Tris-HCl buffer for 30 minutes to expose NHS esters

 Maleimide-C3-NTA
Catalog #: M035

• 0.006 Molarity (M) Maleimide-C3-NTA / Tris-HCl buffer

 Trizma® hydrochloride solution
by Sigma Aldrich
Catalog #: T2319

🕒 00:30:00 Maleimide-C3-NTA

- 3 After washing, the functionalisation molecule is bound to the exposed NHS ester groups by submerging the cantilever in 100 nM Hyaluronic Acid Binding Protein / Tris-HCl buffer solution (Sample 1 HABP/HA and Sample 4 HABP/HAase) or 1% w/v Bovine Serum Albumin / ultra-pure water (Sample 2 BSA/HA) for 1 hour

🧪 0.000001 Molarity (M) Hyaluronic Acid Binding Protein / Tris-HCl buffer

🧪 1 Mass/Volume Percent Bovine Serum Albumin / ultra-pure water

🧪 Hyaluronic Acid Binding Protein, Bovine Nasal Cartilage, Biotinylated
by [Merck Millipore](#)
Catalog #: 385911

🧪 Trizma® hydrochloride solution
by [Sigma Aldrich](#)
Catalog #: T2319

🧪 MilliQ Water

🧪 Bovine Serum Albumin
by [Sigma Aldrich](#)
Catalog #: A9647

🕒 01:00:00 Functionalisation molecule

- 4 The excess maleimide is quenched with 50 mM 2-Mercaptoethanol / ultra-pure water by submerging the cantilevers for 1 minute

🧪 2-Mercaptoethanol
by [Sigma Aldrich](#)
Catalog #: M3148

🧪 MilliQ Water

🕒 00:01:00 2-mercaptoethanol

🧪 0.05 Molarity (M) 2-Mercaptoethanol / ultra-pure water

⚠ SAFETY INFORMATION

Please refer to SDS for 2-Mercaptoethanol

<https://www.sigmaaldrich.com/catalog/product/sigma/m3148?lang=en&ion=GB>

- 5 After a final wash, the functionalised cantilevers are kept submerged in ultra-pure water until mounting on the AFM holder



MilliQ Water



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited