

# Biolistic Transformation of *Amphidinium*

Robert Nimmo, Ellen Nisbet, Adrian Barbrook, Chris Howe.

## Abstract

A method to transform the chloroplast of *Amphidinium carterae* using biolistics.

[orcid.org/0000-0002-2469-0611](https://orcid.org/0000-0002-2469-0611)

[orcid.org/0000-0002-6975-8640](https://orcid.org/0000-0002-6975-8640)

[orcid.org/0000-0003-4487-196X](https://orcid.org/0000-0003-4487-196X)

**Citation:** Robert Nimmo, Ellen Nisbet, Adrian Barbrook, Chris Howe. Biolistic Transformation of *Amphidinium*. **protocols.io**

[dx.doi.org/10.17504/protocols.io.hnmb5c6](https://dx.doi.org/10.17504/protocols.io.hnmb5c6)

**Published:** 13 Apr 2017

## Protocol

To be completed in advance

### Step 1.

Prepare artificial minicircle (pAmpPSBA-Atz) at 1mg ml<sup>-1</sup>. pAmpPSBA-Atz contains a bacterial origin of replication and an ampicillin resistance marker, so transformation and growth using a suitable *E. coli* strain, followed by a miniprep, is the simplest method to achieve this. pAmpPSBA-Atz also contains a linearised minicircle containing the core region and a modified form of PSBA to confer resistance to the herbicide Atrazine.

### Step 2.

Grow wild type *Amphidinium carterae* in *f/2* medium to early log-phase.

### Step 3.

Prepare *f/2* plates (1.5% agarose in *f/2* medium) and leave to dry for 3-4 hours before storing at 4°C. (This drying step is not strictly necessary, but allows the cells to be transferred to the plates much more quickly, which reduces the subsequent time between plating and shooting).

### Step 4.

Autoclave stopping screens, macrocarriers and 1,550 PSI rupture disks (all supplied by BioRad).

## Cell Preparation

### Step 5.

Count the cells. For each plate to be transformed (including necessary controls), between  $1-5 \times 10^7$  cells are needed. Harvest an appropriate volume for the number of plates to be shot.

### Step 6.

Centrifuge the cells at low g (1,500g) and resuspend the pellet in a very small volume (500 $\mu$ l) of fresh f/2 medium.

### Step 7.

Spot the cells onto the centre of the f/2 plate dropwise. Leave to dry. Cells should form a circle of no more than 2cm diameter to maximise the number of cells in the region bombarded by gold particles. Note that *Amphidinium* does not survive long term on solid plates, so aim to shoot as soon after drying as practical.

## DNA Preparation

### Step 8.

*DNA Precipitation is carried out using Seashell Technology's "DNAdel™ gold carrier" delivery system. 550nm gold particles are used throughout.*

### Step 9.

Sonicate the DNAdel™ gold particles (supplied at 50mg ml<sup>-1</sup>) to dissociate aggregates. 1-2 minutes should normally be sufficient, but the presence of aggregates can be detected by eye.

### Step 10.

Dilute the gold particles in the supplied 'binding buffer'\* to a final concentration of 30mg ml<sup>-1</sup>. Use 0.5mg gold per plate for transformation. (e.g. 30ul gold particles and 20ul binding buffer for 3 plates).

\* Supplied by Seashell Technology with the DNAdel™ gold carrier kit.

### Step 11.

Add artificial minicircle to gold particles at a ratio of 2.5 $\mu$ g DNA per mg gold. Using very concentrated DNA stock keeps volume low and allows for high concentrations of DNA and gold to maximise DNA

precipitation onto the gold particles.

## DNA Preparation

### Step 12.

Vortex the DNA-gold particle mixture briefly.

## DNA Preparation

### Step 13.

Add an equal volume of 'precipitation buffer'<sup>\*</sup> to the DNA-gold particle mix and vortex briefly. Allow to stand for 3 minutes.

<sup>\*</sup> Supplied by Seashell Technology with the DNAdel™ gold carrier kit.

### Step 14.

Centrifuge at 10,000g for 10s to pellet the gold particles.

### Step 15.

Remove the supernatant and add 500µl ice cold 100% ethanol. Vortex briefly.

### Step 16.

Centrifuge at 10,000g for 10s to pellet the gold particles.

### Step 17.

Remove the supernatant. Add 10µl ice cold 100% ethanol for each 0.5mg gold in the preparation. (e.g. add 30µl for the 1.5mg gold prepared for shooting 3 plates).

### Step 18.

Briefly sonicate the solution to minimise aggregation and allow for reliable delivery. Again, 1-2 minutes will normally be sufficient, but sonication should continue until even suspension of the gold particles is achieved.

## Bombardment Preparation

### Step 19.

Wash the macrocarriers in 70% ethanol and allow to dry.

### Step 20.

Wash each rupture disk and macrocarrier in isopropanol and allow to dry.

### Step 21.

Transfer a 10µl aliquot of the prepared gold particles (see above) to the centre of each macrocarrier. Allow to dry.

## Bombardment

### Step 22.

*Bombardment is carried out in a Biorad 'Biolistic PDS-1000/He' device using Biorad stopping screens, macrocarriers and rupture disks.*

### Step 23.

Shooting is carried out using 1,550 PSI rupture disks. (These show higher rates of transformation than 1,100 or 1,350 PSI disks.)

### Step 24.

Shooting is carried out in vacuum at  $\geq 25$  in. Hg.

### Step 25.

Plates for transformation are placed in the middle slot in the Biolistics device (*i.e. the highest possible position for a plate. Higher slots are occupied by the rupture disk and macrocarrier*).

## Post Shooting

### Step 26.

*Cells should be resuspended as soon as possible after shooting to minimise cell death.*

### Step 27.

Resuspend all cells from each plate in 5ml fresh *f/2* medium. As many cells as possible should be washed off the surface of the plate and collected in the desired vessel for growth.

### Step 28.

Leave to recover overnight.

### Step 29.

Dilute the cells as and apply selection criteria. (*e.g. each plate's cells are diluted to a final volume of 50ml with 2µg ml<sup>-1</sup> Atrazine*).

### Step 30.

Monitor growth by microscopy every 1-2 days. Untransformed *Amphidinium carterae* are typically dead at 10-12 days post-application of Atrazine.