



Oct 01, 2019

## MojoSort™ Human CD4 T Cell Selection Kit Column Protocol - CD4 Nanobeads [↗](#)

Sam Li<sup>1</sup><sup>1</sup>BioLegend
[1](#) *Works for me* [dx.doi.org/10.17504/protocols.io.7tyhnpw](https://doi.org/10.17504/protocols.io.7tyhnpw)

BioLegend

Sam Li  
BioLegend

### ABSTRACT

BioLegend MojoSort™ nanobeads work in commonly used separation columns, based on our internal research as well as validation by external testing by academic labs. This simple protocol consists of following the MojoSort™ protocol to label the cells with **pre-diluted** MojoSort™ reagents and using the columns as indicated by the manufacturer.

**Note:** Due to the properties of our beads, it may be possible to use far fewer beads than with other commercial suppliers. We recommend a titration to find the best dilution factor. However, as a general rule, dilutions ranging from 1:3 to 1:20 for the Nanobeads can be used. Please contact BioLegend Technical Service ([tech@biolegend.com](mailto:tech@biolegend.com)) if further assistance is needed.

### EXTERNAL LINK

<https://www.biolegend.com/protocols/mojosort-human-cd4-t-cell-selection-kit-column-protocol-cd4-nanobeads/4762/>

### GUIDELINES

MojoSort™ magnetic particles can be used with other commercially available magnetic separators, both free standing magnets and column-based systems. Because MojoSort™ protocols are optimized for the MojoSort™ separator, the protocols may need to be adjusted for other systems. Please contact BioLegend Technical Service ([tech@biolegend.com](mailto:tech@biolegend.com)) for more information and guidance. We do not recommend using MojoSort™ particles for BD's IMag™ or Life Technologies' DynaMag™.

**This protocol is to be used with the two-step Human CD4 T Cell Selection Kit (Cat. No. 480038). First use the CD14 Selection Kit for monocyte depletion, followed by the CD4 Nanobeads to select CD4<sup>+</sup> T Cells.**

### MATERIALS

NAME <a href="#">▼</a>	CATALOG # <a href="#">▼</a>	VENDOR <a href="#">▼</a>
MojoSort™ Buffer	480017	BioLegend
MojoSort™ Human CD4 T Cell Selection Kit	480038	BioLegend

### MATERIALS TEXT

Additional reagents:

- commercially available cell separation columns
- 5 mL polypropylene tubes

- 1 After monocyte depletion using the CD14 selection protocol, centrifuge the pooled unlabeled CD4<sup>+</sup>T Cells at 300xg for 5 minutes, discard supernatant, and resuspend in an appropriate volume of MojoSort™ Buffer. Count and adjust the cell concentration to 1 x 10<sup>8</sup> cells/mL.

**Note:** Keep MojoSort™ Buffer on ice throughout the procedure.

- 2 Aliquot 100 µL (10<sup>7</sup> cells), or desired amount of cells, into a new tube.

- 3 Vortex the antibody-conjugated Nanobeads (to resuspend) at max speed, 5 touches, and prepare the dilutions to test. **Add 10  $\mu\text{L}$  of pre-diluted conjugated Nanobeads.** Mix well and **incubate on ice for 15 minutes**. Scale up the volume accordingly if separating more cells. For example, add 100  $\mu\text{L}$  of pre-diluted Nanobeads for separating  $1 \times 10^8$  cells in 1 ml of MojoSort™ Buffer. When working with less than  $10^7$  cells, use indicated volumes for  $10^7$  cells.
- 4 Add the appropriate amount of MojoSort™ Buffer and proceed to separation. At least 500  $\mu\text{L}$  is needed for column separation.

**Note:** There are several types of commercially available columns, depending on your application. Choose the one that fits best your experiment:

	Max. number of labeled cells	Max. number of total cells	Cell suspension volume	Column rinse volume	Cell wash volume	Elution volume
<b>Small Capacity</b>	$1 \times 10^7$	$2 \times 10^8$	500 $\mu\text{L}$ for up to $10^8$ cells	1ml	1 ml	1 ml
<b>Medium Capacity</b>	$1 \times 10^8$	$2 \times 10^9$	500 $\mu\text{L}$ for up to $10^9$ cells	3ml	3 ml	5 ml
<b>Large Capacity</b>	$1 \times 10^9$	$2 \times 10^{10}$	500 $\mu\text{L}$ for up to $10^{10}$ cells	20-50ml	30 ml	20 ml

Example of magnetic separation with medium capacity columns:

- 5 Place the column in a magnetic separator that fits the column.
- 6 Rinse the column with 3 mL of cell separation buffer.
- 7 Add the labeled cell suspension to the column through a 30  $\mu\text{m}$  filter and collect the fraction containing the unlabeled cells.
- 8 Wash the cells in the column **3 times** with 3 mL of buffer and collect the fraction containing the unlabeled cells. Combine with the collected fraction from step 3. These cells may be useful as controls, to monitor purity/yield, or other purposes.
- 9 Take away the column from the magnet and place it on a tube. Then add 5 mL of buffer and flush out the magnetically labeled fraction with a plunger or supplied device. These are the positively isolated cells of interest; do not discard. To increase the purity of the magnetically labeled fraction repeat the isolation process with a new, freshly prepared column.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited