

1M Tris-HCl

Sarah Hessen-Schmidt

Abstract

Citation: Sarah Hessen-Schmidt 1M Tris-HCl. **protocols.io**

[dx.doi.org/10.17504/protocols.io.fjhbki6](https://doi.org/10.17504/protocols.io.fjhbki6)

Published: 10 Aug 2016


Protocol

Step 1.

Dissolve 30.275g Trizma base in 200ml MilliQ



REAGENTS

 Trizma Base 93362 by Contributed by users

Step 2.

Adjust to pH 7.5 with HCl (35ml)

Step 3.

Adjust total volume to 250ml with MilliQ

Step 4.

Autoclave