OPEN ACCESS



Determination of total glutathione (GSH) and oxidized glutathione (GSSG) levels of RBCs infected with Plasmodium falciparum (Pf RBCs)

Hiroko Asahi

Abstract

Citation: Hiroko Asahi Determination of total glutathione (GSH) and oxidized glutathione (GSSG) levels of RBCs infected with Plasmodium falciparum (Pf RBCs). **protocols.io**

dx.doi.org/10.17504/protocols.io.i35cgq6

Published: 23 Jul 2017

Protocol

Step 1.

step 1

Preparations of PfRBCs for total GSH and GSSG estimation:

1. Harvest PfRBCs*/RBCs (5.4 x 10⁸ cells) in cultures and washed with phosphate-buffered saline.

*at 5% parasitemia

- 2. Disrupt the packed PfRBCs/RBCs by three freeze-thaw cycles.
- 3. Add 100 µl 5% (w/v) sulfosalicylic acid solution (SSA) to the disrupted PfRBCs/RBCs, and vortex.
- 4. Collect the supernatants by centrifugation for 5 min at 8,000g at 4°C, and transfer to 500 μ l chilled, ultrapure water* (final volume 600 μ l).

*Wako Pure Chemical Industries, Osaka, Japan, catalogue no. 210-01303

Step 2.

step 2

Quantification of total glutathione (GSH + GSSG) and GSSG by using the total Glutathione Quantification kit or the GSSG/GSH Quantification kit*.

*total Glutathione Quantification kit (catalogue no. T419) and GSSG/GSH Quantification kit (catalogue no. G257)

Dojindo Molecular Technologies, Inc., Rockville, MD, USA

Step 3.

step 3

Quantification of total glutathione (GSH + GSSG) and GSSG using the GSSG/GSH Quantification kit (G257)

- 1. For total GSH estimation, prepare GSH standard solutions in 0.5% SSA (50, 25.0, 12.5, 6.25, 3.13, 1.57, 0 μ M).
- 2. For GSSG estimation, prepare GSSG standard solutions in 0.5% SSA (25.0, 12.5, 6.25, 3.13, 1.57, 0.78, 0 μ M).
- 3. Add 4 μ l masking solution to 200 μ l GSSG standard solutions, and 200 μ l samples prepared by the procedure above (step 1) in a micro tube, and mix well.
- 4. Place 40 μl of GSSG standards, GSH standards, and samples in a 96-well microtiter plate.
- 5. Add 60 µl Buffer solution to each well.
- 6. Incubate the plate at 37 °C for 1 h.
- 7. Add 60 μl of Substrate working solution to each well.
- 8. Add 60 µl of Enzyme/coenzyme working solution to each well, and mix well.
- 9. Incubate the plate at 37 °C for 10 min.
- 10. Read the absorbance at 405 nm (OD_{405}) on a microplate reader, and extrapolate concentrations of total GSH and GSSG from each standard curve.