



May 08,
2019

Working

U Cinn - Body Composition & Carcass Analysis [↗](#)

Patrick Tso¹, Dana Lee¹

¹University of Cincinnati

[dx.doi.org/10.17504/protocols.io.2qigdue](https://doi.org/10.17504/protocols.io.2qigdue)

Mouse Metabolic Phenotyping Centers

Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

Summary:

Total body composition in live, un-anaesthetized small animals and carcasses will reveal absolute amounts of body fat, lean tissue and body water via a quantitative magnetic resonance (QMR) instrument, EchoMRI, (Echo Medical Systems, LLC, Houston, TX). This instrument uses the differences in the nuclear magnetic resonance properties of hydrogen atoms in organic and non-organic environments to fractionate signals originating from fat, lean tissue and free water.

EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=193&docType=Protocol>

MATERIALS

NAME	CATALOG #	VENDOR	CAS NUMBER	RRID
EchoMRI-100 Whole Body Composition Analyzer for Mice	EchoMRI-100	Echo Medical Systems		
Mouse Restraint	H100-30	Echo Medical Systems		

MATERIALS TEXT

Note:

[EchoMRI](#), [RRID:SCR_017104](#)

- 1 Insert the calibration tube into opening on right side of the EchoMRI-100 as far in as possible.
- 2 Select "**Calibrate**" at the bottom of the screen to calibrate the system.
- 3 After calibration has passed, weigh the animal and carefully place in the restrainer tube.
- 4 Insert the restrainer tube into the opening on the right side of the EchoMRI-100 and:
 - a. Select "**New Experiment**" at the bottom of the screen
 - b. Enter data for the **Group**,
 - c. Enter data for the **Subject**
 - d. **Notes** (Body weight should be included in the "**Notes**" field)
- 5 Select "**Start Experiment**" to start measuring the body composition.

Each run will take approximately 1 minute.

It is recommended that each animal is measured 2 or 3 times to determine the average of the repeated runs.

- 6 When the small box in the upper, left-hand corner reads "Experiment Complete," remove the restrainer from the machine, and return the animal to its home cage.
- 7 Repeat steps 4-7 for all additional animals.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited